



DR4 siRNA (h): sc-35218

BACKGROUND

Tumor necrosis factor (TNF) is a pleiotropic cytokine whose function is mediated by two distinct cell surface receptors, designated TNF-R1 and TNF-R2, which are expressed on most cell types. TNF function is primarily mediated through TNF-R1 signaling. Both receptors belong to the growing TNF receptor superfamily which includes FAS antigen and CD40. TNF-R1 contains a cytoplasmic motif, termed the "death domain", that has been found to be necessary for the transduction of the apoptotic signal. The death domain is also found in several other receptors, including FAS, DR2 (or TRUNDD), DR3 (death receptor 3), DR4 and DR5. TRUNDD, DR4 and DR5 are receptors for the apoptosis-inducing cytokine TRAIL. A non-death domain-containing receptor, designated decoy receptor (DcR1 or TRID), also specifically associates with TRAIL and may play a role in cellular resistance to apoptotic stimuli.

CHROMOSOMAL LOCATION

Genetic locus: TNFRSF10A (human) mapping to 8p21.3.

PRODUCT

DR4 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see DR4 shRNA Plasmid (h): sc-35218-SH and DR4 shRNA (h) Lentiviral Particles: sc-35218-V as alternate gene silencing products.

For independent verification of DR4 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-35218A, sc-35218B and sc-35218C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

DR4 siRNA (h) is recommended for the inhibition of DR4 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

DR4 (B-9): sc-8411 is recommended as a control antibody for monitoring of DR4 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor DR4 gene expression knockdown using RT-PCR Primer: DR4 (h)-PR: sc-35218-PR (20 μ l, 497 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

- Mendoza, F.J., et al. 2008. Death receptor-4 (DR4) expression is regulated by transcription factor NF κ B in response to etoposide treatment. *Apoptosis* 13: 756-770.
- Seon, M.R., 2010. Isoangustone A present in hexane/ethanol extract of *Glycyrrhiza uralensis* induces apoptosis in DU145 human prostate cancer cells via the activation of DR4 and intrinsic apoptosis pathway. *Mol. Nutr. Food Res.* 54: 1329-1339.
- Martínez-Paniagua, M.A., et al. 2011. Mcl-1 and YY1 inhibition and induction of DR5 by the BH3-mimetic obatoclax (GX15-070) contribute in the sensitization of B-NHL cells to TRAIL apoptosis. *Cell Cycle* 10: 2792-2805.
- Jo, M., et al. 2012. Anti-cancer effect of bee venom toxin and melittin in ovarian cancer cells through induction of death receptors and inhibition of JAK2/Stat3 pathway. *Toxicol. Appl. Pharmacol.* 258: 72-81.
- Choi, K.E., et al. 2014. Myriocin induces apoptotic lung cancer cell death via activation of DR4 pathway. *Arch. Pharm. Res.* 37: 501-511.
- Hou, L., et al. 2015. Effect of vitamin E succinate on the expression of the tumor necrosis factor-related apoptosis-inducing ligand (TRAIL) receptor in gastric cancer cells and CD4⁺ T cells. *Mol. Biosyst.* 11: 3119-3128.
- Lim, S.C., et al. 2016. Involvement of DR4/JNK pathway-mediated autophagy in acquired TRAIL resistance in Hep G2 cells. *Int. J. Oncol.* 49: 1983-1990.
- Kim, B.R., et al. 2019. RUNX3 enhances TRAIL-induced apoptosis by upregulating DR5 in colorectal cancer. *Oncogene* 38: 3903-3918.
- Yeh, L.Y., et al. 2020. The miR-372-ZBTB7A oncogenic axis suppresses TRAIL-R2 associated drug sensitivity in oral carcinoma. *Front. Oncol.* 10: 47.
- Zhang, Y.Y., et al. 2022. Licochalcone B induces DNA damage, cell cycle arrest, apoptosis, and enhances TRAIL sensitivity in hepatocellular carcinoma cells. *Chem. Biol. Interact.* 365: 110076.
- Li, L., et al. 2022. Decitabine enhances the tumoricidal potential of TRAIL via the epigenetic regulation of death receptor 4 in gastric cancer. *J. Gastrointest. Oncol.* 13: 2799-2808.

RESEARCH USE

For research use only, not for use in diagnostic procedures.