

GADD 153 siRNA (m): sc-35438

BACKGROUND

GADD 153 has been described as a growth arrest and DNA damage-inducible gene that encodes a C/EBP-related nuclear protein. This protein has also been designated C/EBP-homologous protein (CHOP-10). GADD 153 expression is induced by a variety of cellular stresses, inducing nutrient deprivation and metabolic perturbations. GADD 153 functions to block cells in G₁ to S phase in cell cycle progression and acts by dimerizing with other C/EBP proteins to direct GADD 153 dimers away from "classical" C/EBP binding sites, recognizing instead unique "nonclassical" sites. Thus GADD 153 acts as a negative modulator of C/EBP-like proteins in certain terminally differentiated cells, similar to the regulatory function of Id on the activity of MyoD and MyoD-related proteins involved in the development of muscle cells.

CHROMOSOMAL LOCATION

Genetic locus: Ddit3 (mouse) mapping to 10 D3.

PRODUCT

GADD 153 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GADD 153 shRNA Plasmid (m): sc-35438-SH and GADD 153 shRNA (m) Lentiviral Particles: sc-35438-V as alternate gene silencing products.

For independent verification of GADD 153 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-35438A, sc-35438B and sc-35438C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GADD 153 siRNA (m) is recommended for the inhibition of GADD 153 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

GADD 153 (B-3): sc-7351 is recommended as a control antibody for monitoring of GADD 153 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GADD 153 gene expression knockdown using RT-PCR Primer: GADD 153 (m)-PR: sc-35438-PR (20 μ l, 443 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

- Li, K., et al. 2013. Arsenic trioxide alleviates airway hyperresponsiveness and promotes apoptosis of CD4⁺ T lymphocytes: evidence for involvement of the ER stress-CHOP pathway. *Ir. J. Med. Sci.* 182: 573-583.
- Yuan, Y., et al. 2015. The roles of oxidative stress, endoplasmic reticulum stress, and autophagy in aldosterone/mineralocorticoid receptor-induced podocyte injury. *Lab. Invest.* 95: 1374-1386.
- Marwarha, G., et al. 2016. Palmitate-induced endoplasmic reticulum stress and subsequent C/EBP α homologous protein activation attenuates leptin and Insulin-like growth factor 1 expression in the brain. *Cell. Signal.* 28: 1789-1805.
- Rao, Z., et al. 2017. Hyperglycemia aggravates hepatic ischemia and reperfusion injury by inhibiting liver-resident macrophage M2 polarization via C/EBP homologous protein-mediated endoplasmic reticulum stress. *Front. Immunol.* 8: 1299.
- Kim, W.H., et al. 2019. C/EBP homologous protein deficiency inhibits statin-induced myotoxicity. *Biochem. Biophys. Res. Commun.* 508: 857-863.
- Fan, C., et al. 2020. Melatonin suppresses ER stress-dependent proapoptotic effects via AMPK in bone mesenchymal stem cells during mitochondrial oxidative damage. *Stem Cell Res. Ther.* 11: 442.
- Suzuki, R., et al. 2020. Intracellular accumulation of advanced glycation end-products induces osteoblast apoptosis via endoplasmic reticulum stress. *J. Bone Miner. Res.* 35: 1992-2003.
- Yang, Z., et al. 2020. Heme induces BECN1/ATG5-mediated autophagic cell death via ER stress in neurons. *Neurotox. Res.* 38: 1037-1048.
- Molagoda, I.M.N., et al. 2021. Anthocyanins from *Hibiscus syriacus* L. inhibit NLRP3 inflammasome in BV2 microglia cells by alleviating NF κ B- and ER stress-induced Ca²⁺ accumulation and mitochondrial Ros production. *Oxid. Med. Cell. Longev.* 2021: 1246491.
- Nazir, L.A., et al. 2022. Synthesis and anti-melanoma effect of 3-O-prenyl glycyrrhetic acid against B16F10 cells via induction of endoplasmic reticulum stress-mediated autophagy through ERK/AKT signaling pathway. *Front. Oncol.* 12: 890299.

RESEARCH USE

For research use only, not for use in diagnostic procedures.