

CLOSTRIDIUM PERFRINGENS TYPES C & D-TETANUS TOXOID

Colorado Serum
Toxoid

Clostridium Perfringens Types C & D-Tetanus Toxoid

Contents: Purified formalin detoxified filtrates of highly toxic cultures of Clostridium perfringens types C and D micro-organisms and of Tetanus Toxin.

Contains thimerosal as a preservative.

Indications: For the vaccination of healthy cattle, sheep, and goats as an aid in the prevention of enterotoxemia caused by Clostridium perfringens types C and D. For the vaccination of healthy swine as an aid in the prevention of enterotoxemia caused by Clostridium perfringens type C. For the long-term protection against tetanus.

Dosage and Administration: Inject 2 mL subcutaneously or intramuscularly. Repeat the full dose in 3 to 4 weeks.

Precaution(s): Shake well. Each dose must have a proportionate share of the precipitate for proper responses. Use the entire contents when first opened.

Store in dark at 2° to 7°C. Do not freeze.

CAUTION(s): Anaphylactoid reaction sometimes follows administration of products of this nature. If noted, administer adrenaline or equivalent.

For veterinary use only.

WARNING(s): Do not vaccinate with 21 days before slaughter.

Professional Biological Company
4950 York St., P.O. Box 16428
Denver, CO 80216-0428
1-800-525-2065
www.professionalbiological.com

Disclaimer

Every effort has been made by Jeffers to ensure the accuracy of the information listed above. However, it remains your responsibility to become familiar with the products you are purchasing. Please consult your family veterinarian if you have any questions.