

Bo-Bac™ 2X Antiserum

Boehringer Ingelheim

Antiserum

Actinomyces Pyogenes-Escherichia Coli-Pasteurella Multocida-Salmonella Typhimurium Antibody,
Bovine Origin

Composition: The concentrated serum is prepared from the blood of cattle hyperimmunized with Actinomyces pyogenes, Pasteurella multocida Carter's Serotype A, Escherichia coli serotype 78:K80:NM, and Salmonella typhimurium.

Contains 0.2% cresol as a preservative.

Indications: Recommended for the prophylaxis and treatment of diseases caused by Actinomyces pyogenes, Pasteurella multocida, Escherichia coli and Salmonella typhimurium in cattle.

Dosage and Administration: Prophylactic: 15 mL per 50 lbs. of body weight administered subcutaneously as soon as possible after birth.

Therapeutic: 30 mL per 50 lbs. of body weight, repeat every 12 to 24 hours (depending upon condition of the animal) until improvement is satisfactory.

Precaution(s): Store out of direct sunlight at a temperature between 35–45°F (2–7°C). Avoid freezing. Shake well before using. Use the entire contents when first opened.

Caution(s): Anaphylactoid reactions may occur.

Antidote(s): Epinephrine.

Warning(s): Do not vaccinate within 21 days before slaughter.

Do not vaccinate with bovine virus diarrhea vaccine, bovine rhinotracheitis vaccine or bovine parainfluenza3 vaccine, or Haemophilus somnus or Pasteurella haemolytica bacterins within 21 days after use of this serum.

Boehringer Ingelheim Vetmedica, Inc.
2621 North Belt Highway
St. Joseph, MO 64506–2002
1–800–325–9167
www.bi-vetmedica.com

Disclaimer

Every effort has been made by Jeffers to ensure the accuracy of the information listed above. However, it remains your responsibility to become familiar with the products you are purchasing. Please consult your family veterinarian if you have any questions.