



p70 S6 kinase α siRNA (h): sc-36165

BACKGROUND

In studies to elucidate key regulatory pathways in signal transduction, several protein serine/threonine (Ser/Thr) kinases have been identified, including two distinct families of 40S ribosomal protein S6 Ser/Thr kinases present in somatic animal cells, designated p70 S6 kinase and p90 Rsk kinase. p90 Rsk kinase is maximally activated within minutes of addition of growth factors or phorbol ester to cultured cells followed by activation of p70 S6 kinase. Both enzymes are regulated by serine/threonine phosphorylation, suggesting that specific kinases may exist upstream in the signaling pathway that regulate these kinases. In fact, evidence suggests that one such family of activating enzymes includes the members of the ERK MAP kinase family. The ERK MAP kinases are, in turn, regulated by phosphorylation at threonine and tyrosine residues by a protein kinase designated MEK.

CHROMOSOMAL LOCATION

Genetic locus: RPS6KB1 (human) mapping to 17q23.1.

PRODUCT

p70 S6 kinase α siRNA (h) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see p70 S6 kinase α shRNA Plasmid (h): sc-36165-SH and p70 S6 kinase α shRNA (h) Lentiviral Particles: sc-36165-V as alternate gene silencing products.

For independent verification of p70 S6 kinase α (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-36165A and sc-36165B.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

p70 S6 kinase α siRNA (h) is recommended for the inhibition of p70 S6 kinase α expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

p70 S6 kinase α (H-9): sc-8418 is recommended as a control antibody for monitoring of p70 S6 kinase α gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor p70 S6 kinase α gene expression knockdown using RT-PCR Primer: p70 S6 kinase α (h)-PR: sc-36165-PR (20 μ l, 533 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

- Vazquez-Martin, A., et al. 2009. The antidiabetic drug metformin suppresses HER2 (erbB-2) oncoprotein overexpression via inhibition of the mTOR effector p70S6K1 in human breast carcinoma cells. *Cell Cycle* 8: 88-96.
- Hisamatsu, T., et al. 2013. Potential role of mTORC2 as a therapeutic target in clear cell carcinoma of the ovary. *Mol. Cancer Ther.* 12: 1367-1377.
- Eckerd, F., et al. 2014. Regulatory effects of a Mnk2-eIF4E feedback loop during mTORC1 targeting of human medulloblastoma cells. *Oncotarget* 5: 8442-8451.
- Chatterjee, A., et al. 2015. Rapamycin-induced G₁ cell cycle arrest employs both TGF- β and Rb pathways. *Cancer Lett.* 360: 134-140.
- Hac, A., et al. 2015. S6K1 controls autophagosome maturation in autophagy induced by sulforaphane or serum deprivation. *Eur. J. Cell Biol.* 94: 470-481.
- Yun, S.M., et al. 2016. Melatonin enhances arsenic trioxide-induced cell death via sustained upregulation of Redd1 expression in breast cancer cells. *Mol. Cell. Endocrinol.* 422: 64-73.
- Lu, Z., et al. 2016. Downregulation of p70S6K enhances cell sensitivity to rapamycin in esophageal squamous cell carcinoma. *J. Immunol. Res.* 2016: 7828916.
- Chen, G.H., et al. 2017. Inhibition of miR-128-3p by tongxinluo protects human cardiomyocytes from ischemia/reperfusion injury via upregulation of p70s6k1/p-p70s6k1. *Front. Pharmacol.* 8: 775.
- Chaturvedi, N.K., et al. 2018. Improved therapy for medulloblastoma: targeting hedgehog and PI3K-mTOR signaling pathways in combination with chemotherapy. *Oncotarget* 9: 16619-16633.
- Gao, Q., et al. 2019. Distinct role of 4E-BP1 and S6K1 in regulating autophagy and hepatitis B virus (HBV) replication. *Life Sci.* 220: 1-7.
- Chen, G., et al. 2020. Cardiomyocyte-derived small extracellular vesicles can signal eNOS activation in cardiac microvascular endothelial cells to protect against ischemia/reperfusion injury. *Theranostics* 10: 11754-11774.

RESEARCH USE

For research use only, not for use in diagnostic procedures.