

VR1 siRNA (m): sc-36827

BACKGROUND

Vanilloid receptor 1 (VR1), also designated capsaicin receptor, is a nonselective cation channel, structurally related to members of the TRP family of ion channels. VR1 is activated by capsaicin, the active ingredient in chili peppers, by heat and by an increase in protons at sites of infection, inflammation and ischemia. By creating moderately acidic conditions, protons are able to lower the temperature threshold for VR1 activation, thus identifying VR1 as a molecular integrator of chemical and physical stimuli that elicit pain. VR1 is expressed in primary sensory neurons and vagal nerves and activated VR1 induces the influx of cations, particularly Ca^{2+} and Na^{+} ions. The vanilloid receptor may also be a molecular target for endogenous anandamide, in addition to the cannabinoid receptors, in the nervous and cardiovascular systems.

REFERENCES

1. Caterina, M.J., et al. 1997. The capsaicin receptor: a heat-activated ion channel in the pain pathway. *Nature* 389: 816-824.
2. Tominaga, M., et al. 1998. The cloned capsaicin receptor integrates multiple pain-producing stimuli. *Neuron* 21: 531-543.
3. Cesare, P., et al. 1999. Ion channels gated by heat. *PNAS* 96: 7658-7663.
4. Sasamura, T. and Kuraishi, Y. 1999. Peripheral and central actions of capsaicin and VR1 receptor. *Jpn. J. Pharmacol.* 80: 275-280.
5. Zygmunt, P.M., et al. 1999. Vanilloid receptors on sensory nerves mediate the vasodilator action of anandamide. *Nature* 400: 452-457.
6. Jin, Y.H., et al. 2004. Purinergic and vanilloid receptor activation releases glutamate from separate cranial afferent terminals in nucleus tractus solitarius. *J. Neurosci.* 24: 4709-4717.
7. Jung, J., et al. 2004. Phosphorylation of vanilloid receptor 1 by Ca^{2+} /calmodulin-dependent kinase II regulates its vanilloid binding. *J. Biol. Chem.* 279: 7048-7054.
8. Karai, L.J., et al. 2004. Vanilloid receptor 1 regulates multiple calcium compartments and contributes to Ca^{2+} -induced Ca^{2+} release in sensory neurons. *J. Biol. Chem.* 279: 16377-16387.

CHROMOSOMAL LOCATION

Genetic locus: *Trpv1* (mouse) mapping to 11 B4.

PRODUCT

VR1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μM solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see VR1 shRNA Plasmid (m): sc-36827-SH and VR1 shRNA (m) Lentiviral Particles: sc-36827-V as alternate gene silencing products.

For independent verification of VR1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-36827A, sc-36827B and sc-36827C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20°C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20°C , avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μl of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μl of RNase-free water makes a 10 μM solution in a 10 μM Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

VR1 siRNA (m) is recommended for the inhibition of VR1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μM in 66 μl . Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

VR1 (E-8): sc-398417 is recommended as a control antibody for monitoring of VR1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor VR1 gene expression knockdown using RT-PCR Primer: VR1 (m)-PR: sc-36827-PR (20 μl , 542 bp). Annealing temperature for the primers should be $55-60^{\circ}\text{C}$ and the extension temperature should be $68-72^{\circ}\text{C}$.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.