

# MerTK siRNA (h): sc-37127

## BACKGROUND

MerTK, also called c-Mer, is a member of the Mer/Axl/Tyro3 receptor kinase family. It is a 984 residue transmembrane protein made up of one tyrosine kinase domain, two Fibronectin type-III domains and two immunoglobulin-like C2-type domains. MerTK is the mammalian ortholog of the chicken retroviral oncogene product v-Eyk. This protein plays a critical role in macrophage activation, platelet aggregation, clot stability and the efficient removal of apoptotic cells. Specifically, MerTK acts as a signaling molecule, triggering outer segment ingestion in the retinal pigment epithelium (RPE) phagocytic process. Evidence suggests that MerTK signals via interaction with phosphatidylinositol-specific phospholipase C  $\gamma 2$  (PI-PLC  $\gamma 2$ ). When the gene encoding for MerTK is mutated, the RPE phagocytosis pathway is disrupted and autosomal recessive retinitis pigmentosa (RP) may result, leading to degeneration of retinal photoreceptor cells.

## REFERENCES

1. Graham, D.K., et al. 1994. Cloning and mRNA expression analysis of a novel human protooncogene, c-Mer. *Cell Growth Differ.* 5: 647-657.
2. Gal, A., et al. 2000. Mutations in MerTK, the human orthologue of the RCS rat retinal dystrophy gene, cause retinitis pigmentosa. *Nat. Genet.* 26: 270-271.
3. D'Cruz, P.M., et al. 2000. Mutation of the receptor tyrosine kinase gene MerTK in the retinal dystrophic RCS rat. *Hum. Mol. Genet.* 9: 645-651.
4. Kumar, A., et al. 2001. Retinitis pigmentosa: mutations in a receptor tyrosine kinase gene, MerTK. *J. Biosci.* 26: 3-5.

## CHROMOSOMAL LOCATION

Genetic locus: MERTK (human) mapping to 2q13.

## PRODUCT

MerTK siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see MerTK shRNA Plasmid (h): sc-37127-SH and MerTK shRNA (h) Lentiviral Particles: sc-37127-V as alternate gene silencing products.

For independent verification of MerTK (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-37127A, sc-37127B and sc-37127C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

MerTK siRNA (h) is recommended for the inhibition of MerTK expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

MerTK (B-1): sc-365499 is recommended as a control antibody for monitoring of MerTK gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended:

1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor MerTK gene expression knockdown using RT-PCR Primer: MerTK (h)-PR: sc-37127-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## SELECT PRODUCT CITATIONS

1. Roggia, M.F. and Ueta, T. 2015.  $\alpha_v\beta_5$  Integrin/FAK/PGC-1 $\alpha$  pathway confers protective effects on retinal pigment epithelium. *PLoS ONE* 10: e0134870.
2. Jhang, K.A., et al. 2018. Sulforaphane rescues amyloid- $\beta$  peptide-mediated decrease in MerTK expression through its anti-inflammatory effect in human THP-1 macrophages. *J. Neuroinflammation* 15: 75.
3. Li, Y., et al. 2019. The role of endothelial MERTK during the inflammatory response in lungs. *PLoS ONE* 14: e0225051.
4. Ge, X., et al. 2022. Chromium (VI)-induced ALDH1A1/EGF axis promotes lung cancer progression. *Clin. Transl. Med.* 12: e1136.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.