



# PGI2 synthase siRNA (m): sc-37237

## BACKGROUND

Prostacyclin (also known as prostaglandin I<sub>2</sub>) is a potent vasodilator and inhibitor of platelet aggregation. The enzyme PGI<sub>2</sub> synthase (also known as prostacyclin synthase) catalyzes the isomerization of prostaglandin H<sub>2</sub> (PGH<sub>2</sub>) to prostacyclin. Although it has absorbance spectral features characteristic of the cytochrome P450s, PGIS has no monooxygenase activity and does not require an external source of electrons to initiate its enzyme reaction. PGI<sub>2</sub> synthase is the single member of family 8 of the cytochrome P450 superfamily. PGI<sub>2</sub> synthase is a polypeptide of 500 amino acids with sequence homology to cholesterol 7- $\alpha$ -monooxygenase, a member of the CYP7 family of cytochrome P450s. The gene which encodes PGI<sub>2</sub> synthase maps to human chromosome 20q13.13.

## REFERENCES

1. Miyata, A., Hara, S., Yokoyama, C., Inoue, H., Ullrich, V. and Tanabe, T. 1994. Molecular cloning and expression of human prostacyclin synthase. *Biochem. Biophys. Res. Commun.* 200: 1728-1734.
2. Wang, L.H. and Chen, L. 1996. Organization of the gene encoding human prostacyclin synthase. *Biochem. Biophys. Res. Commun.* 226: 631-637.
3. Nelson, D.R., Koymans, L., Kamataki, T., Stegeman, J.J., Feyereisen, R., Waxman, D.J., Waterman, M.R., Gotoh, O., Coon, M.J., Estabrook, R.W., Gunsalus, I.C. and Nebert, D.W. 1996. P450 superfamily: update on new sequences, gene mapping, accession numbers and nomenclature. *Pharmacogenetics* 6: 1-42.
4. Yokoyama, C., Yabuki, T., Inoue, H., Tone, Y., Hara, S., Hatae, T., Nagata, M., Takahashi, E.I. and Tanabe, T. 1996. Human gene encoding prostacyclin synthase (PTGIS): genomic organization, chromosomal localization, and promoter activity. *Genomics* 36: 296-304.
5. LocusLink Report (LocusID: 601699). <http://www.ncbi.nlm.nih.gov/LocusLink/>

## CHROMOSOMAL LOCATION

Genetic locus: Ptgis (mouse) mapping to 2 H3.

## PRODUCT

PGI<sub>2</sub> synthase siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see PGI<sub>2</sub> synthase shRNA Plasmid (m): sc-37237-SH and PGI<sub>2</sub> synthase shRNA (m) Lentiviral Particles: sc-37237-V as alternate gene silencing products.

For independent verification of PGI<sub>2</sub> synthase (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-37237A, sc-37237B and sc-37237C.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

PGI<sub>2</sub> synthase siRNA (m) is recommended for the inhibition of PGI<sub>2</sub> synthase expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor PGI<sub>2</sub> synthase gene expression knockdown using RT-PCR Primer: PGI<sub>2</sub> synthase (m)-PR: sc-37237-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## SELECT PRODUCT CITATIONS

1. Barbieri, S.S., Amadio, P., Gianellini, S., Zacchi, E., Weksler, B.B. and Tremoli, E. 2011. Tobacco smoke regulates the expression and activity of microsomal prostaglandin E synthase-1: role of prostacyclin and NADPH-oxidase. *FASEB J.* 25: 3731-3740.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.