

Smad2/3 siRNA (h): sc-37238

BACKGROUND

Smad proteins, the mammalian homologs of the *Drosophila* mothers against decapentaplegic (Mad), have been implicated as downstream effectors of TGF β /BMP signaling. Smad1 (also designated Madr1 or JV4-1) and Smad5 are effectors of BMP-2 and BMP-4 function, while Smad2 (also designated Madr2 or JV18-1) and Smad3 are involved in TGF β and Activin-mediated growth modulation. Smad4 (also designated DPC4) has been shown to mediate all of the above activities through interaction with various Smad family members. Smad6 and Smad7 regulate the response to Activin/TGF β signaling by interfering with TGF β -mediated phosphorylation of other Smad proteins.

CHROMOSOMAL LOCATION

Genetic locus: SMAD2 (human) mapping to 18q21.1, SMAD3 (human) mapping to 15q22.33.

PRODUCT

Smad2/3 siRNA (h) is a pool of 4 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Smad2/3 shRNA Plasmid (h): sc-37238-SH and Smad2/3 shRNA (h) Lentiviral Particles: sc-37238-V as alternate gene silencing products.

For independent verification of Smad2/3 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-37238A, sc-37238B, sc-37238C and sc-37238D.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Smad2/3 siRNA (h) is recommended for the inhibition of Smad2/3 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Smad2/3 (C-8): sc-133098 is recommended as a control antibody for monitoring of Smad2/3 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

SELECT PRODUCT CITATIONS

1. Zhao, L., et al. 2010. LPS-induced epithelial-mesenchymal transition of intrahepatic biliary epithelial cells. *J. Surg. Res.* 171: 819-825.
2. Han, L. and Gotlieb, A.I. 2011. Fibroblast growth factor-2 promotes *in vitro* mitral valve interstitial cell repair through transforming growth factor- β /Smad signaling. *Am. J. Pathol.* 178: 119-127.
3. Zhao, W.J. and Wang, K. 2013. NGX6 expression improves the sensitivity of tamoxifen-resistant MCF7 cells through modulation of the Smad signaling pathway. *Int. J. Oncol.* 42: 2060-2068.
4. Kim, M.R., et al. 2015. TGF β 1 protects cells from γ -IR by enhancing the activity of the NHEJ repair pathway. *Mol. Cancer Res.* 13: 319-329.
5. Zhang, Y., et al. 2015. Apigenin induces dermal collagen synthesis via Smad2/3 signaling pathway. *Eur. J. Histochem.* 59: 2467.
6. Lee, J., et al. 2016. TGF- β 1 accelerates the DNA damage response in epithelial cells via Smad signaling. *Biochem. Biophys. Res. Commun.* 476: 420-425.
7. Da, C., et al. 2016. Nobiletin inhibits epithelial-mesenchymal transition of human non-small cell lung cancer cells by antagonizing the TGF- β 1/Smad3 signaling pathway. *Oncol. Rep.* 35: 2767-2774.
8. Liu, F., et al. 2018. Metformin ameliorates TGF- β 1-induced osteoblastic differentiation of human aortic valve interstitial cells by inhibiting β -catenin signaling. *Biochem. Biophys. Res. Commun.* 500: 710-716.
9. Wa, Q., et al. 2018. Downregulation of miR-19a-3p promotes invasion, migration and bone metastasis via activating TGF- β signaling in prostate cancer. *Oncol. Rep.* 39: 81-90.
10. Jung, M.Y., et al. 2018. Fatty Acid Synthase is required for profibrotic TGF- β signaling. *FASEB J.* 32: 3803-3815.
11. Fenizia, C., et al. 2018. SMYD3 promotes the epithelial-mesenchymal transition in breast cancer. *Nucleic Acids Res.* 47: 1278-1293.
12. Jung, H., et al. 2018. Fucoxanthin inhibits myofibroblast differentiation and extracellular matrix production in nasal polyp-derived fibroblasts via modulation of Smad-dependent and Smad-independent signaling pathways. *Mar. Drugs* 16: 323.
13. Park, J.S., et al. 2019. Targeting of dermal myofibroblasts through death receptor 5 arrests fibrosis in mouse models of scleroderma. *Nat. Commun.* 10: 1128.

RESEARCH USE

For research use only, not for use in diagnostic procedures.