

CIDE-B siRNA (h): sc-37441

BACKGROUND

The DNA fragmentation factor (DFF) is involved in the caspase-3 apoptotic pathway. DFF is composed of two subunits: DFF-45, also designated ICAD (inhibitor of CAD); and CPAN (caspase-activated nuclease), also designated CAD (caspase-activated deoxyribonuclease). CPAN is a DNase that is responsible for DNA degradation during apoptosis. CPAN is inhibited by DFF-45. Caspase-3 acts to dissociate CPAN from DFF-45, allowing CPAN to enter the nucleus and degrade DNA. CIDE-A and CIDE-B have been identified as proteins that share homology with the N-terminal region of DFF-45. Like CPAN, CIDE-A and CIDE-B promote cell death and DNA fragmentation and are inhibited by DFF-45.

REFERENCES

1. Inohara, N., et al. 1998. CIDE, a novel family of cell death activators with homology to the 45 kDa subunit of the DNA fragmentation factor. *EMBO J.* 17: 2526-2533.
2. Inohara, N., et al. 1999. Identification of regulatory and catalytic domains in the apoptosis nuclease DFF40/CAD. *J. Biol. Chem.* 274: 270-274.
3. Erdtmann, L., et al. 2003. The hepatitis C virus NS2 protein is an inhibitor of CIDE-B-induced apoptosis. *J. Biol. Chem.* 278: 18256-18264.
4. Machado, J.G., et al. 2005. Gene expression profiling of jejunal Peyer's patches in juvenile and adult pigs. *Mamm. Genome* 16: 599-612.
5. Turpaev, K., et al. 2005. Analysis of differentially expressed genes in nitric oxide-exposed human monocytic cells. *Free Radic. Biol. Med.* 38: 1392-1400.
6. Novikova, S.I., et al. 2005. Cocaine-induced changes in the expression of apoptosis-related genes in the fetal mouse cerebral wall. *Neurotoxicol. Teratol.* 27:3-14.
7. Da, L., et al. 2006. Dual promoters control the cell-specific expression of the human cell death-inducing DFF45-like effector B gene. *Biochem. J.* 393: 779-788.
8. Wang, J.M., et al. 2006. Cell death-inducing DFF45-like effector may take part in neuronal apoptosis of the lumbar spinal cord after sciatic nerve injury caused by a firearm. *Mil. Med.* 171: 793-799.

CHROMOSOMAL LOCATION

Genetic locus: CIDEB (human) mapping to 14q12.

PRODUCT

CIDE-B siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CIDE-B shRNA Plasmid (h): sc-37441-SH and CIDE-B shRNA (h) Lentiviral Particles: sc-37441-V as alternate gene silencing products.

For independent verification of CIDE-B (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-37441A, sc-37441B and sc-37441C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

CIDE-B siRNA (h) is recommended for the inhibition of CIDE-B expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

CIDE-B (5X): sc-101244 is recommended as a control antibody for monitoring of CIDE-B gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CIDE-B gene expression knockdown using RT-PCR Primer: CIDE-B (h)-PR: sc-37441-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.