

Hus1 siRNA (m): sc-37546

BACKGROUND

DNA damage or incomplete replication of DNA results in inhibition of cell cycle progression at the G₁-S or G₂-M checkpoints by conserved regulatory mechanisms. Chk1, Rad9 and Hus1 are involved in regulation of cell cycle arrest at the G₂ checkpoint. Chk1 functions as an essential component in the G₂ DNA damage checkpoint by phosphorylating Cdc25C in response to DNA damage, which inhibits mitosis. Hus1 and Rad9 exhibit conserved function in fission yeast and higher eukaryotes. Hus1 has been shown to be phosphorylated in response to DNA damage, a process which requires rad checkpoint genes. Rad9 is thought to be a candidate tumor suppressor gene because it is localized to human chromosome 11q13.2, which is a region containing a number of tumor suppressor loci.

REFERENCES

1. Carr, A.M., et al. 1995. The Chk1 pathway is required to prevent mitosis following cell-cycle arrest at "start". *Curr. Biol.* 5: 1179-1190.
2. Lieberman, H.B., et al. 1996. A human homolog of the *Schizosaccharomyces pombe* Rad9⁺ checkpoint control gene. *Proc. Natl. Acad. Sci. USA* 93: 13890-13895.
3. Sanchez, Y., et al. 1997. Conservation of the Chk1 checkpoint pathway in mammals: linkage of DNA damage to Cdk regulation through Cdc25. *Science* 277: 1497-1501.
4. Peng, C.Y., et al. 1997. Mitotic and G₂ checkpoint control: regulation of 14-3-3 protein binding by phosphorylation of Cdc25C on serine-216. *Science* 277: 1501-1505.
5. O'Connell, M.J., et al. 1997. Chk1 is a Wee1 kinase in the G₂ DNA damage checkpoint inhibiting Cdc2 by Y15 phosphorylation. *EMBO J.* 16: 545-554.
6. Kostrub, C.F., et al. 1998. Hus1p, a conserved fission yeast checkpoint protein, interacts with Rad1p and is phosphorylated in response to DNA damage. *EMBO J.* 17: 2055-2066.
7. Wang, X., et al. 2004. Involvement of Hus1 in the chain elongation step of DNA replication after exposure to camptothecin or ionizing radiation. *Nucleic Acids Res.* 32: 767-775.

CHROMOSOMAL LOCATION

Genetic locus: Hus1 (mouse) mapping to 11 A1.

PRODUCT

Hus1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Hus1 shRNA Plasmid (m): sc-37546-SH and Hus1 shRNA (m) Lentiviral Particles: sc-37546-V as alternate gene silencing products.

For independent verification of Hus1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-37546A, sc-37546B and sc-37546C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Hus1 siRNA (m) is recommended for the inhibition of Hus1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Hus1 (G-3): sc-166440 is recommended as a control antibody for monitoring of Hus1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Hus1 gene expression knockdown using RT-PCR Primer: Hus1 (m)-PR: sc-37546-PR (20 μ l, 520 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.