

CARP siRNA (m): sc-37732

BACKGROUND

A proposed genetic marker of cardiac hypertrophy, CARP (cardiac ankyrin repeat protein) is a nuclear protein with an established role in regulation of cardiac gene expression. A distinct increase in CARP expression occurs in rats with abdominal aorta constriction, spontaneous hypertension and Dahl salt-sensitivity. In cardiomyocytes, CARP inhibits transcription of both cardiac troponin C and atrial natriuretic factor. Specifically, expression of the CARP gene, which lies downstream of the cardiac homeobox gene *Nkx2.5*, inhibits *Nkx2.5* transactivation of atrial natriuretic factor promoter. An increase in CARP expression is observed in the ventricular tissue of patients with end-stage heart failure. The major Ca^{2+} binding protein of cardiac sarcoplasmic reticulum (SR), Calsequestrin (CSQ), upregulates the CARP gene, and may contribute to the development of cardiac hypertrophy and fibrosis. TGF β induces CARP expression in vascular smooth muscle cells (VSMCs), wherein CARP may mediate the inhibitory effects of TGF β on VSMC proliferation.

REFERENCES

- Jeyaseelan, R., et al. 1997. A novel cardiac-restricted target for doxorubicin. CARP, a nuclear modulator of gene expression in cardiac progenitor cells and cardiomyocytes. *J. Biol. Chem.* 272: 22800-22808.
- Zou, Y., et al. 1997. CARP, a cardiac ankyrin repeat protein, is downstream in the *Nkx2-5* homeobox gene pathway. *Development* 124: 793-804.
- Aihara, Y., et al. 2000. Cardiac ankyrin repeat protein is a novel marker of cardiac hypertrophy: role of M-CAT element within the promoter. *Hypertension* 36: 48-53.
- Bang, M.L., et al. 2001. Myopalladin, a novel 145-kilodalton sarcomeric protein with multiple roles in Z-disc and I-band protein assemblies. *J. Cell Biol.* 153: 413-427.
- Kanai, H., et al. 2001. Transforming growth factor- β /Smads signaling induces transcription of the cell type-restricted ankyrin repeat protein CARP gene through CAGA motif in vascular smooth muscle cells. *Circ. Res.* 88: 30-36.
- Zolk, O., et al. 2002. Cardiac ankyrin repeat protein, a negative regulator of cardiac gene expression, is augmented in human heart failure. *Biochem. Biophys. Res. Commun.* 293: 1377-1382.

CHROMOSOMAL LOCATION

Genetic locus: *Ankrd1* (mouse) mapping to 19 C2.

PRODUCT

CARP siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CARP shRNA Plasmid (m): sc-37732-SH and CARP shRNA (m) Lentiviral Particles: sc-37732-V as alternate gene silencing products.

For independent verification of CARP (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-37732A, sc-37732B and sc-37732C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

CARP siRNA (m) is recommended for the inhibition of CARP expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

CARP (G-2): sc-365056 is recommended as a control antibody for monitoring of CARP gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz MarkerTM Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CARP gene expression knockdown using RT-PCR Primer: CARP (m)-PR: sc-37732-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.