

Pax-2 siRNA (h): sc-38745

BACKGROUND

Pax genes contain paired domains with strong homology to genes in *Drosophila* which are involved in programming early development. The PAX2 gene is expressed in primitive cells of the kidney, ureter, eye, ear, and central nervous system. More specifically, in human embryo sections, PAX2 is expressed in the optic vesicle and later in the retina, in the otic vesicle and later in the semicircular canals of the inner ear, and in mesonephros, metanephros, adrenals, spinal cord, and hindbrain. PAX2 mutations can be responsible for renal hypoplasia, either isolated or associated with various ophthalmologic manifestations ranging from retinal coloboma to microphthalmia. Lesions in the PAX6 gene accounts for most cases of aniridia, a congenital malformation of the eye, chiefly characterized by iris hypoplasia, which can cause blindness. PAX6 is involved in other anterior segment malformations besides aniridia, such as Peters anomaly, a major error in the embryonic development of the eye with corneal clouding with variable iridolenticulocorneal adhesions.

REFERENCES

1. Ferrell, R.E., et al. 1980. Autosomal dominant aniridia: probable linkage to acid phosphatase-1 locus on chromosome 2. *Proc. Natl. Acad. Sci. USA* 77: 1580-1582.
2. Hanson, I.M., et al. 1993. PAX6 mutations in aniridia. *Hum. Mol. Genet.* 2: 915-920.
3. Hanson, I.M., et al. 1994. Mutations at the PAX6 locus are found in heterogeneous anterior segment malformations including Peters' anomaly. *Nat. Genet.* 6: 168-173.
4. Narahara, K., et al. 1997. Localization of a 10q breakpoint within the PAX2 gene in a patient with a *de novo* t(10;13) translocation and optic nerve coloboma-renal disease. *J. Med. Genet.* 34: 213-216.

CHROMOSOMAL LOCATION

Genetic locus: PAX2 (human) mapping to 10q24.31.

PRODUCT

Pax-2 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Pax-2 shRNA Plasmid (h): sc-38745-SH and Pax-2 shRNA (h) Lentiviral Particles: sc-38745-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Pax-2 siRNA (h) is recommended for the inhibition of Pax-2 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Pax-2 (60-P): sc-130387 is recommended as a control antibody for monitoring of Pax-2 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Pax-2 gene expression knockdown using RT-PCR Primer: Pax-2 (h)-PR: sc-38745-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

1. Browne, J.A., et al. 2014. Open chromatin mapping identifies transcriptional networks regulating human epididymis epithelial function. *Mol. Hum. Reprod.* 20: 1198-1207.
2. Liu, P., et al. 2015. Upregulation of PAX2 promotes the metastasis of esophageal cancer through interleukin-5. *Cell. Physiol. Biochem.* 35: 740-754.
3. Browne, J.A., et al. 2018. Region-specific innate antiviral responses of the human epididymis. *Mol. Cell. Endocrinol.* 473: 72-78.

RESEARCH USE

For research use only, not for use in diagnostic procedures.