

VSX1 siRNA (m): sc-38807

BACKGROUND

Like other "paired-like" homeodomain family members, the visual system homeobox gene 1 (VSX1) is instrumental in craniofacial and ocular development; VSX1 plays a distinct role in retinal development. Also known as RINX (retinal inner nuclear layer homeobox), the VSX1 gene is expressed in embryonic craniofacial structures and in the adult retina. VSX1 is abundantly expressed in the inner nuclear layer (INL) of the retina. In mice, *Vsx1* first detected in the bipolar cells of the retina five days postnatal. The VSX1 gene is also expressed in WERI, a retinoblastoma cell line that expresses retinal cone genes. The human VSX1 gene maps to chromosome 20p11.21 and encodes a 365 amino acid protein with five known splice variants. VSX1 mutations are implicated in two distinct corneal dystrophies, posterior polymorphous dystrophy (PPD) and keratoconus.

REFERENCES

1. Semina, E.V., Mintz-Hittner, H.A. and Murray, J.C. 2000. Isolation and characterization of a novel human paired-like homeodomain-containing transcription factor gene, VSX1, expressed in ocular tissues. *Genomics* 63: 289-293.
2. Hayashi, T., Huang, J. and Deeb, S.S. 2000. RINX (VSX1), a novel homeobox gene expressed in the inner nuclear layer of the adult retina. *Genomics* 67: 128-139.
3. Chow, R.L., Snow, B., Novak, J., Looser, J., Freund, C., Vidgen, D., Ploder, L. and McInnes, R.R. 2001. VSX1, a rapidly evolving paired-like homeobox gene expressed in cone bipolar cells. *Mech. Dev.* 109: 315-322.
4. Heon, E., Greenberg, A., Kopp, K.K., Rootman, D., Vincent, A.L., Billingsley, G., Priston, M., Dorval, K.M., Chow, R.L., McInnes, R.R., Heathcote, G., Westall, C., Sutphin, J.E., Semina, E., Bremner, R. and Stone, E.M. 2002. VSX1: a gene for posterior polymorphous dystrophy and keratoconus. *Hum. Mol. Genet.* 11: 1029-1036.
5. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 605020. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: *Vsx1* (mouse) mapping to 2 G3.

PRODUCT

VSX1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see VSX1 shRNA Plasmid (m): sc-38807-SH and VSX1 shRNA (m) Lentiviral Particles: sc-38807-V as alternate gene silencing products.

For independent verification of VSX1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-38807A, sc-38807B and sc-38807C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

VSX1 siRNA (m) is recommended for the inhibition of VSX1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

VSX1 (G-11): sc-393699 is recommended as a control antibody for monitoring of VSX1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor VSX1 gene expression knockdown using RT-PCR Primer: VSX1 (m)-PR: sc-38807-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.