

FXR siRNA (h): sc-38848

BACKGROUND

The steroid receptor superfamily acts through direct association with DNA sequences known as hormone response elements (HREs) and binds DNA as either homo- or heterodimers. The promiscuous mediator of heterodimerization, RXR, is the receptor for 9-*cis* retinoic acid, and dimerizes with VDR, TR, PPAR, and several novel receptors including LXR (also referred to as RLD-1) and FXR. FXR and LXR fall into a category of proteins termed "orphan receptors" because of their lack of a defined function, and in the case of LXR, the lack of a defined ligand. FXR has been shown to bind a class of lipid molecules called farnesoids. LXR/RXR heterodimers have highest affinity for DR-4 DNA elements while FXR/RXR heterodimers bind IR-1 elements. Both LXR/RXR and FXR/RXR heterodimers retain their responsiveness to 9-*cis* retinoic acid.

CHROMOSOMAL LOCATION

Genetic locus: NR1H4 (human) mapping to 12q23.1.

PRODUCT

FXR siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see FXR shRNA Plasmid (h): sc-38848-SH and FXR shRNA (h) Lentiviral Particles: sc-38848-V as alternate gene silencing products.

For independent verification of FXR (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-38848A, sc-38848B and sc-38848C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

FXR siRNA (h) is recommended for the inhibition of FXR expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

FXR (D-3): sc-25309 is recommended as a control antibody for monitoring of FXR gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor FXR gene expression knockdown using RT-PCR Primer: FXR (h)-PR: sc-38848-PR (20 μ l, 599 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

- Wang, Y.D., et al. 2015. Farnesoid X receptor antagonizes JNK signaling pathway in liver carcinogenesis by activating SOD3. *Mol. Endocrinol.* 29: 322-331.
- Joshi, S., et al. 2016. Bile acids-mediated overexpression of MUC4 via FAK-dependent c-Jun activation in pancreatic cancer. *Mol. Oncol.* 10: 1063-1077.
- Xu, W., et al. 2017. Dihydroartemisinin protects against alcoholic liver injury through alleviating hepatocyte steatosis in a farnesoid X receptor-dependent manner. *Toxicol. Appl. Pharmacol.* 315: 23-34.
- Liu, Z.Q., et al. 2018. Ciliogenesis is reciprocally regulated by PPARA and NR1H4/FXR through controlling autophagy *in vitro* and *in vivo*. *Autophagy* 14: 1011-1027.
- Qiao, P., et al. 2018. Farnesoid X receptor inhibits proliferation of human colorectal cancer cells via the miR-135A1/CCNG2 signaling pathway. *Oncol. Rep.* 40: 2067-2078.
- Wu, Z.Y., et al. 2019. Farnesoid X receptor agonist GW4064 indirectly inhibits HCV entry into cells via down-regulating scavenger receptor class B type I. *Eur. J. Pharmacol.* 853: 111-120.
- Mao, J., et al. 2020. Effects and mechanism of the bile acid (farnesoid X) receptor on the Wnt/ β -catenin signaling pathway in colon cancer. *Oncol. Lett.* 20: 337-345.
- Lee, Y.J., et al. 2020. The role of nuclear receptor subfamily 1 group H member 4 (NR1H4) in colon cancer cell survival through the regulation of c-Myc stability. *Mol. Cells* 43: 459-468.
- Chaudhari, S.N., et al. 2021. A microbial metabolite remodels the gut-liver axis following bariatric surgery. *Cell Host Microbe* 29: 408-424.e7.
- Fu, Y., et al. 2022. Alisol B 23-acetate adjusts bile acid metabolism via hepatic FXR-BSEP signaling activation to alleviate atherosclerosis. *Phytomedicine* 101: 154120.
- Xu, X.W., et al. 2022. Selective photoaffinity probe for monitoring farnesoid X receptor expression in cultured cells. *Anal. Chem.* 94: 10722-10729.

RESEARCH USE

For research use only, not for use in diagnostic procedures.