## SANTA CRUZ BIOTECHNOLOGY, INC.

# hnRNP F/H (G-8): sc-390085



### BACKGROUND

Heterogeneous nuclear ribonucleoproteins (hnRNPs) constitute a set of polypeptides that contribute to pre-mRNA processing and transport, and also bind heterogeneous nuclear RNA (hnRNA), which are the transcripts produced by RNA polymerase II. hnRNP complexes are the major constituents of the spliceosome and in particular, the hnRNP A1 protein is one of the major pre-mRNA/mRNA binding proteins and also one of the most abundant proteins in the nucleus. hnRNP A1 and A2/B1 regulate the processing of premRNA by directly antagonizing the association of various splicing factors and by influencing the splice site selection on pre-mRNA. The majority of hnRNP proteins components are localized to the nucleus; however some shuttle between the nucleus and the cytoplasm. Most hnRNP proteins, including hnRNP C1 and C2, contain one or more RNA binding domains and are implicated in the processing of pre-mRNA. hnRNPs F and H are largely related factors that preferentially associate with poly(rG) regions on RNA. Isoforms of these proteins are often generated by alternative processing of the premRNA and by posttranslational modifications such as phosphorylation on serines and threonines and methylation of arginines.

#### REFERENCES

- 1. Swanson, M.S., et al. 1987. Primary structure of human nuclear ribonucleoprotein particle C proteins. Mol. Cell. Biol. 7: 1731-1739.
- 2. Gorlach, M., et al. 1994. The determinants of RNA-binding specificity of the heterogeneous nuclear ribonucleoprotein C proteins. J. Biol. Chem. 269: 23074-23078.
- 3. Honore, B., et al. 1995. Heterogeneous nuclear ribonucleoproteins H, H', and F are members of a ubiquitously expressed subfamily of related but distinct proteins encoded by genes mapping to different chromosomes. J. Biol. Chem. 270: 28780-28789.
- 4. Badolato, J., et al. 1995. Identification and characterization of a novel human RNA-binding protein. Gene 166: 323-327.
- 5. Siomi, H., et al. 1995. A nuclear localization domain in the hnRNP A1 protein. J. Cell Biol. 129: 551-560.
- 6. Hanamura, A., et al. 1998. Regulated tissue-specific expression of antagonistic pre-mRNA splicing factors. RNA 4: 430-444.
- 7. Melcak, I., et al. 2000. Nuclear pre-mRNA compartmentalization: trafficking of released transcripts to splicing factor reservoirs. Mol. Biol. Cell 11: 497-510.

## SOURCE

hnRNP F/H (G-8) is a mouse monoclonal antibody raised against amino acids 1-300 of hnRNP F/H of human origin.

## **PRODUCT**

Each vial contains 200  $\mu$ g lgG<sub>1</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## **RESEARCH USE**

For research use only, not for use in diagnostic procedures.

#### **APPLICATIONS**

hnRNP F/H (G-8) is recommended for detection of hnRNP F, hnRNP H1 and hnRNP H' of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Molecular Weight of hnRNP F/H: 48 kDa.

Positive Controls: K-562 nuclear extract: sc-2130 or Hep G2 cell lysate: sc-2227 or Jurkat whole cell lysate: sc-2204.

#### **RECOMMENDED SUPPORT REAGENTS**

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-lgG K BP-HRP: sc-516102 or m-lgG K BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgGk BP-FITC: sc-516140 or m-IgGk BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850. 4) Immunohistochemistry: use m-lgG $\kappa$  BP-HRP: sc-516102 with DAB, 50X: sc-24982 and Immunohistomount: sc-45086, or Organo/Limonene Mount: sc-45087.

### DATA





hnRNP F/H (G-8): sc-390085. Western blot analysis of hnBNP F/H expression in K-562 nuclear extract (A) and Hep G2 (B), Jurkat (C), A-375 (D) and DU 145 (E) whole cell lysates.

hnRNP F/H (G-8): sc-390085. Immunoperoxidase staining of formalin fixed, paraffin-embedded human duodenum tissue showing nuclear staining of glandular cells

## SELECT PRODUCT CITATIONS

- 1. Russo, A., et al. 2010. hnRNP H1 and intronic G runs in the splicing control of the human rpL3 gene. Biochim. Biophys. Acta 1799: 419-428.
- 2. Russo, A., et al. 2011. Autoregulatory circuit of human rpL3 expression requires hnRNP H1, NPM and KHSRP. Nucleic Acids Res. 39: 7576-7585.

#### **STORAGE**

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.