

# TTC12 (B-1): sc-390229

## BACKGROUND

The tetratricopeptide repeat (TPR) motif is a degenerate, 34 amino acid sequence found in many proteins and acts to mediate protein-protein interactions in various pathways. At the sequence level, there can be up to 16 tandem TPR repeats, each of which has a helix-turn-helix shape that stacks on other TPR repeats to achieve ligand binding specificity. TTC12 (tetratricopeptide repeat domain 12), also known as TPARM, is a 705 amino acid cytoplasmic protein containing a TPR repeat domain and three armadillo repeat motifs. TTC12 is expressed in testis, prostate, lung, germinal center B cells, neuroblastoma, teratocarcinoma, colon cancer and gastric cancer. The gene encoding TTC12 is located in a region of human chromosome 11 that is commonly deleted in a variety of cancers, therefore, it is considered a candidate tumor suppressor gene. Human chromosome 11 houses over 1,400 genes and comprises nearly 4% of the human genome.

## REFERENCES

1. Katoh, M., et al. 2003. Identification and characterization of TPARM gene in silico. *Int. J. Oncol.* 23: 1213-1217.
2. Gelernter, J., et al. 2006. Haplotype spanning TTC12 and ANKK1, flanked by the DRD2 and NCAM1 loci, is strongly associated to nicotine dependence in two distinct American populations. *Hum. Mol. Genet.* 15: 3498-3507.

## CHROMOSOMAL LOCATION

Genetic locus: TTC12 (human) mapping to 11q23.2; Ttc12 (mouse) mapping to 9 A5.3.

## SOURCE

TTC12 (B-1) is a mouse monoclonal antibody specific for an epitope mapping between amino acids 7-41 near the N-terminus of TTC12 of human origin.

## PRODUCT

Each vial contains 200 µg IgG<sub>1</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

TTC12 (B-1) is available conjugated to agarose (sc-390229 AC), 500 µg/0.25 ml agarose in 1 ml, for IP; to HRP (sc-390229 HRP), 200 µg/ml, for WB, IHC(P) and ELISA; to either phycoerythrin (sc-390229 PE), fluorescein (sc-390229 FITC), Alexa Fluor® 488 (sc-390229 AF488), Alexa Fluor® 546 (sc-390229 AF546), Alexa Fluor® 594 (sc-390229 AF594) or Alexa Fluor® 647 (sc-390229 AF647), 200 µg/ml, for WB (RGB), IF, IHC(P) and FCM; and to either Alexa Fluor® 680 (sc-390229 AF680) or Alexa Fluor® 790 (sc-390229 AF790), 200 µg/ml, for Near-Infrared (NIR) WB, IF and FCM.

Blocking peptide available for competition studies, sc-390229 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% stabilizer protein).

Alexa Fluor® is a trademark of Molecular Probes, Inc., Oregon, USA

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## APPLICATIONS

TTC12 (B-1) is recommended for detection of TTC12 of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

TTC12 (B-1) is also recommended for detection of TTC12 in additional species, including bovine and porcine.

Suitable for use as control antibody for TTC12 siRNA (h): sc-96505, TTC12 siRNA (m): sc-154750, TTC12 shRNA Plasmid (h): sc-96505-SH, TTC12 shRNA Plasmid (m): sc-154750-SH, TTC12 shRNA (h) Lentiviral Particles: sc-96505-V and TTC12 shRNA (m) Lentiviral Particles: sc-154750-V.


Molecular Weight of TTC12: 79 kDa.

Positive Controls: TTC12 (h): 293T Lysate: sc-114915, human testis extract: sc-363781 or rat testis extract: sc-2400.


## RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## DATA



TTC12 (B-1): sc-390229. Western blot analysis of TTC12 expression in non-transfected: sc-117752 (A) and human TTC12 transfected: sc-114915 (B) 293T whole cell lysates and human testis (C) and rat testis (D) tissue extracts.



TTC12 (B-1): sc-390229. Immunofluorescence staining of methanol-fixed HeLa cells showing cytoplasmic localization.

## SELECT PRODUCT CITATIONS

1. Chen, W., et al. 2022. Biallelic mutations of TTC12 and TTC21B were identified in Chinese patients with multisystem ciliopathy syndromes. *Hum. Genomics* 16: 48.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.