

ETO-2 (G-7): sc-390332



The Power to Question

BACKGROUND

ETO and ETO-2, which are alternatively designated MTG8 and MTG16, respectively, are members of the ETO transcription factor family. These transcription factors are characterized by a zinc-finger domain and four conserved domains, of which domain II is required for dimerization between family members. ETO and ETO-2 may function to mediate interactions between DNA binding proteins and transcriptional regulators, such as N-CoR. Frequently, the t(8;21) translocation of ETO produces the AML-1/ETO oncoprotein, which consists of the first 177 amino acids of AML-1 and all but the first 30 amino acids of ETO. AML-1/ETO expression is observed in 12-15% of acute myelogenous, M2 subtype leukemias. The AML-1/ETO fusion proteins associate with multimeric N-CoR/mSin3/HDAC1 complexes, block differentiation and induce transcriptional repression by altering chromatin remodeling.

REFERENCES

1. Erickson, P.F., et al. 1994. The ETO portion of acute myeloid leukemia t(8;21) fusion transcript encodes a highly evolutionarily conserved, putative transcription factor. *Cancer Res.* 54: 1782-1786.
2. Erickson, P.F., et al. 1996. ETO and AML1 phosphoproteins are expressed in CD34⁺ hematopoietic progenitors: implications for t(8;21) leukemogenesis and monitoring residual disease. *Blood* 88: 1813-1823.
3. Wolford, J.K., et al. 1998. Structure and expression of the human MTG8/ETO gene. *Gene* 212: 103-109.
4. Wang, J., et al. 1998. ETO, fusion partner in t(8;21) acute myeloid leukemia, represses transcription by interaction with the human N-CoR/mSin3/HDAC1 complex. *Proc. Natl. Acad. Sci. USA* 95: 10860-10865.
5. Westendorf, J.J., et al. 1998. The t(8;21) fusion product, AML-1-ETO, associates with C/EBP- α , inhibits C/EBP- α -dependent transcription, and blocks granulocytic differentiation. *Mol. Cell. Biol.* 18: 322-333.
6. Davis, J.N., et al. 1999. ETO-2, a new member of the ETO-family of nuclear proteins. *Oncogene* 18: 1375-1383.

CHROMOSOMAL LOCATION

Genetic locus: CBFA2T3 (human) mapping to 16q24.3; Cbfa2t3 (mouse) mapping to 8 E1.

SOURCE

ETO-2 (G-7) is a mouse monoclonal antibody specific for an epitope mapping between amino acids 586-605 near the C-terminus of ETO-2 of mouse origin.

PRODUCT

Each vial contains 200 μ g IgM kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin. Also available as TransCruz reagent for Gel Supershift and ChIP applications, sc-390332 X, 200 μ g/0.1 ml.

Blocking peptide available for competition studies, sc-390332 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% stabilizer protein).

APPLICATIONS

ETO-2 (G-7) is recommended for detection of ETO-2 of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for ETO-2 siRNA (h): sc-35344, ETO-2 siRNA (m): sc-35345, ETO-2 shRNA Plasmid (h): sc-35344-SH, ETO-2 shRNA Plasmid (m): sc-35345-SH, ETO-2 shRNA (h) Lentiviral Particles: sc-35344-V and ETO-2 shRNA (m) Lentiviral Particles: sc-35345-V.

ETO-2 (G-7) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.


Molecular Weight of ETO-2: 76 kDa.

Positive Controls: ETO-2 (h): 293T Lysate: sc-116819 or SUP-T1 whole cell lysate: sc-364796.


RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein L-Agarose: sc-2336 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

DATA



ETO-2 (G-7): sc-390332. Western blot analysis of ETO-2 expression in non-transfected: sc-117752 (A) and human ETO-2 transfected: sc-116819 (B) 293T whole cell lysates.



ETO-2 (G-7): sc-390332. Western blot analysis of ETO-2 expression in SUP-T1 whole cell lysate.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.