

# DDX33 (F-5): sc-390436

## BACKGROUND

DEAD-box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp, are putative RNA helicases implicated in several cellular processes involving modifications of RNA secondary structure and ribosome/spliceosome assembly. Based on their distribution patterns, some members of this family may be involved in embryogenesis, spermatogenesis and cellular growth and division. DDX33 (DEAD (Asp-Glu-Ala-His) box polypeptide 33), also known as DHX33, is a 707 amino acid nucleolar protein belonging to the DEAD box helicase family. Containing a helicase ATP-binding domain and a helicase C-terminal domain, DDX33 is encoded by a gene located on human chromosome 17. Chromosome 17 comprises over 2.5% of the human genome and encodes over 1,200 genes. Two isoforms of DDX33 exists due to alternative splicing events.

## REFERENCES

- Schmid, S.R. and Linder, P. 1992. D-E-A-D protein family of putative RNA helicases. *Mol. Microbiol.* 6: 283-291.
- Will, C.L., Urlaub, H., Achsel, T., Gentzel, M., Wilm, M. and Lührmann, R. 2002. Characterization of novel SF3b and 17S U2 snRNP proteins, including a human Prp5p homologue and an SF3b DEAD-box protein. *EMBO J.* 21: 4978-4988.
- Abdelhaleem, M., Maltais, L. and Wain, H. 2003. The human DDX and DHX gene families of putative RNA helicases. *Genomics* 81: 618-622.
- Xu, Y.Z., Newnham, C.M., Kameoka, S., Huang, T., Konarska, M.M. and Query, C.C. 2004. Prp5 bridges U1 and U2 snRNPs and enables stable U2 snRNP association with intron RNA. *EMBO J.* 23: 376-385.
- Cordin, O., Tanner, N.K., Doère, M., Linder, P. and Banroques, J. 2004. The newly discovered Q motif of DEAD-box RNA helicases regulates RNA-binding and helicase activity. *EMBO J.* 23: 2478-2487.
- Andersen, J.S., Lam, Y.W., Leung, A.K., Ong, S.E., Lyon, C.E., Lamond, A.I. and Mann, M. 2005. Nucleolar proteome dynamics. *Nature* 433: 77-83.
- Linder, P. 2006. DEAD-box proteins: a family affair—active and passive players in RNP-remodeling. *Nucleic Acids Res.* 34: 4168-4180.

## CHROMOSOMAL LOCATION

Genetic locus: DHX33 (human) mapping to 17p13.2; Dhx33 (mouse) mapping to 11 B4.

## SOURCE

DDX33 (F-5) is a mouse monoclonal antibody specific for an epitope mapping between amino acids 621-655 near the C-terminus of DDX33 of human origin.

## PRODUCT

Each vial contains 200 µg IgG<sub>2b</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-390436 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% stabilizer protein).

## APPLICATIONS

DDX33 (F-5) is recommended for detection of DDX33 of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

DDX33 (F-5) is also recommended for detection of DDX33 in additional species, including equine.

Suitable for use as control antibody for DDX33 siRNA (h): sc-93804, DDX33 siRNA (m): sc-143037, DDX33 shRNA Plasmid (h): sc-93804-SH, DDX33 shRNA Plasmid (m): sc-143037-SH, DDX33 shRNA (h) Lentiviral Particles: sc-93804-V and DDX33 shRNA (m) Lentiviral Particles: sc-143037-V.

Molecular Weight (predicted) of DDX33 isoforms: 79/60 kDa.

Molecular Weight (observed) of DDX33 isoforms: 69/80 kDa.


Positive Controls: 3T3-L1 cell lysate: sc-2243, HeLa nuclear extract: sc-2120 or HeLa whole cell lysate: sc-2200.

## RECOMMENDED SUPPORT REAGENTS


To ensure optimal results, the following support reagents are recommended:

- Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.
- Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).
- Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## DATA



DDX33 (F-5): sc-390436. Western blot analysis of DDX33 expression in HeLa whole cell lysate (A) and HeLa nuclear extract (B).



DDX33 (F-5): sc-390436. Western blot analysis of DDX33 expression in 3T3-L1 whole cell lysate.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.