

CALM (G-8): sc-390580

BACKGROUND

Clathrin-coated pits and vesicles are assembled for receptor-mediated endocytosis through interaction with clathrin associated protein complexes. Vesicle transport is mediated from the *trans*-Golgi network by the adapter complex AP-1 and from the plasma membrane by the AP-2 complex. The AP-1 and AP-2 adapter protein complexes consist of clathrin binding adaptin proteins (γ and $\beta 1$ for AP-1, α and $\beta 2$ for AP-2) and two smaller subunits known as AP50 and AP17. The α - and β -adaptin chains have a similar two-domain organization with C-terminal domains that vary in both sequence and length. α -adaptin splice variants A and C display variable relative expression levels and differential distribution in different tissues. AP180 (also designated AP-3 or F1-20) is a synapse-specific clathrin assembly protein. The protein CALM (clathrin assembly protein lymphoid myeloid leukemia) is highly homologous to AP180 and may also be involved in clathrin assembly.

REFERENCES

1. Robinson, M.S. 1989. Cloning of cDNAs encoding two related 100 kD coated vesicle proteins (α -adaptins). *J. Cell Biol.* 108: 833-842.
2. Kirchhausen, T., et al. 1989. Structural and functional division into two domains of the large (100 to 115 kDa) chains of the clathrin-associated protein complex AP-2. *Proc. Natl. Acad. Sci. USA* 86: 2612-2616.
3. Robinson, M.S. 1990. Cloning and expression of γ -adaptin, a component of Clathrin-coated vesicles associated with the Golgi apparatus. *J. Cell Biol.* 111: 2319-2326.
4. Ponnambalam, S., et al. 1990. Conservation and diversity in families of coated vesicle adaptins. *J. Biol. Chem.* 265: 4814-4820.
5. Morris, S.A., et al. 1993. Clathrin assembly protein AP180: primary structure, domain organization and identification of a clathrin binding site. *EMBO J.* 12: 667-675.
6. Ball, C.L., et al. 1995. Expression and localization of α -adaptin isoforms. *J. Cell Sci.* 108: 2865-2875.
7. Mellman, I. 1996. Endocytosis and molecular sorting. *Annu. Rev. Cell Dev. Biol.* 12: 575-625.

CHROMOSOMAL LOCATION

Genetic locus: PICALM (human) mapping to 11q14.2; Picalm (mouse) mapping to 7 E1.

SOURCE

CALM (G-8) is a mouse monoclonal antibody raised against amino acids 287-420 mapping within an internal region of CALM of human origin.

PRODUCT

Each vial contains 200 μ g IgG₁ kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

CALM (G-8) is recommended for detection of CALM of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

CALM (G-8) is also recommended for detection of CALM in additional species, including canine.

Suitable for use as control antibody for CALM siRNA (h): sc-29882, CALM siRNA (m): sc-29883, CALM shRNA Plasmid (h): sc-29882-SH, CALM shRNA Plasmid (m): sc-29883-SH, CALM shRNA (h) Lentiviral Particles: sc-29882-V and CALM shRNA (m) Lentiviral Particles: sc-29883-V.


Molecular Weight of CALM: 62-72 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200, NIH/3T3 whole cell lysate: sc-2210 or K-562 whole cell lysate: sc-2203.


RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850. 4) Immunohistochemistry: use m-IgG κ BP-HRP: sc-516102 with DAB, 50X: sc-24982 and Immunohistomount: sc-45086, or Organo/Limonene Mount: sc-45087.

DATA



CALM (G-8): sc-390580. Western blot analysis of CALM expression in HeLa (A), K-562 (B), NIH/3T3 (C), 3T3-L1 (D) and RAW 264.7 (E) whole cell lysates.



CALM (G-8): sc-390580. Immunoperoxidase staining of formalin fixed, paraffin-embedded human small intestine tissue showing cytoplasmic and membrane staining of glandular cells.

RESEARCH USE

For research use only, not for use in diagnostic procedures.