

## PP2A-A $\beta$ siRNA (m): sc-39180

### BACKGROUND

The protein phosphatase (PP) holoenzyme is a trimeric complex composed of a regulatory subunit, a variable subunit and a catalytic subunit. Four major families of protein phosphatase catalytic subunits have been identified, designated PP1, PP2A, PP2B (calcineurin) and PP2C. An additional protein phosphatase catalytic subunit, PPX (also known as PP4) is a putative member of a novel PP family. The PP2A family comprises subfamily members PP2A $\alpha$  and PP2A $\beta$ . The PP2A catalytic subunit associates with a variety of regulatory subunits. Regulatory subunits include PP2A-A $\alpha$  and -A $\beta$ , PP2A-B $\alpha$  and -B $\beta$ , PP2A-C $\alpha$  and -C $\beta$ , and PP2A-B56 $\alpha$  and -B56 $\beta$ . PP2A-A $\beta$ , also known as PPP2R1B, is a 601 amino acid protein that contains 15 HEAT repeats. Each HEAT repeat consists of two  $\alpha$  helices joined by a hydrophilic region, known as the intrarepeat loop. The repeat units may be arranged laterally to form a rod-like structure.

### REFERENCES

- Hemmings, B.A., et al. 1990.  $\alpha$ - and  $\beta$ -forms of the 65-kDa subunit of protein phosphatase 2A have a similar 39 amino acid repeating structure. *Biochemistry* 29: 3166-3173.
- Baysal, B.E., et al. 1998. Genomic organization and precise physical location of protein phosphatase 2A regulatory subunit A  $\beta$  isoform gene on chromosome band 11q23. *Gene* 217: 107-116.
- Wang, S.S., et al. 1998. Alterations of the PPP2R1B gene in human lung and colon cancer. *Science* 282: 284-287.
- Campbell, I.G. and Manolitsas, T. 1999. Absence of PPP2R1B gene alterations in primary ovarian cancers. *Oncogene* 18: 6367-6369.
- Takagi, Y., et al. 2000. Alterations of the PPP2R1B gene located at 11q23 in human colorectal cancers. *Gut* 47: 268-271.
- Hemmer, S., et al. 2002. Alterations in the suppressor gene PPP2R1B in parathyroid hyperplasias and adenomas. *Cancer Genet. Cytogenet.* 134: 13-17.
- Lubert, E.J. and Sarge, K.D. 2003. Interaction between protein phosphatase 2A and members of the importin beta superfamily. *Biochem. Biophys. Res. Commun.* 303: 908-913.

### CHROMOSOMAL LOCATION

Genetic locus: Ppp2r1b (mouse) mapping to 9 A5.3.

### PRODUCT

PP2A-A $\beta$  siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see PP2A-A $\beta$  shRNA Plasmid (m): sc-39180-SH and PP2A-A $\beta$  shRNA (m) Lentiviral Particles: sc-39180-V as alternate gene silencing products.

For independent verification of PP2A-A $\beta$  (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-39180A, sc-39180B and sc-39180C.

### STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

### APPLICATIONS

PP2A-A $\beta$  siRNA (m) is recommended for the inhibition of PP2A-A $\beta$  expression in mouse cells.

### SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

### GENE EXPRESSION MONITORING

PP2A-A $\alpha$ / $\beta$  (4G7): sc-13600 is recommended as a control antibody for monitoring of PP2A-A $\beta$  gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

### RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor PP2A-A $\beta$  gene expression knockdown using RT-PCR Primer: PP2A-A $\beta$  (m)-PR: sc-39180-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.