

APIP (G-2): sc-393194

BACKGROUND

APIP (Apaf-1-interacting protein), also known as APIP2, is a member of the aldolase class II family and has a highly conserved C-terminal from *C. elegans* to humans. It is ubiquitously expressed, with high expression levels in heart, kidney and skeletal muscle. Alternative splicing produces two isoforms of APIP. Isoform 1 is the full length, 242 amino acid protein; isoform 2 is missing residues 1-38 and contains a distinct sequence for amino acids 39-53. APIP plays an important role in preventing muscle ischemic damage. It suppresses hypoxia-induced cell death by inducing the activation of Akt and ERK 1/2, which are responsible for inhibition of caspase-9 via phosphorylation, and competing with caspase-9 to bind the caspase recruitment domain (CARD) of Apaf-1. Through these mechanisms, APIP negatively regulates the activation of caspase-9 and Apaf-1-mediated cell death.

REFERENCES

1. Lai, C.H., et al. 2000. Identification of novel human genes evolutionarily conserved in *Caenorhabditis elegans* by comparative proteomics. *Genome Res.* 10: 703-713.
2. Cao, G., et al. 2004. Cloning of a novel Apaf-1-interacting protein: a potent suppressor of apoptosis and ischemic neuronal cell death. *J. Neurosci.* 24: 6189-6201.
3. Cho, D.H., et al. 2004. Induced inhibition of ischemic/hypoxic injury by APIP, a novel Apaf-1-interacting protein. *J. Biol. Chem.* 279: 39942-39950.
4. Anichini, A., et al. 2006. Apaf-1 signaling in human melanoma. *Cancer Lett.* 238: 168-179.
5. Schafer, Z.T. and Kornbluth, S. 2006. The apoptosome: physiological, developmental, and pathological modes of regulation. *Dev. Cell* 10: 549-561.

CHROMOSOMAL LOCATION

Genetic locus: APIP (human) mapping to 11p13; Apip (mouse) mapping to 2 E2.

SOURCE

APIP (G-2) is a mouse monoclonal antibody raised against amino acids 171-242 mapping at the C-terminus of APIP of human origin.

PRODUCT

Each vial contains 200 µg IgG_{2a} kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APIP (G-2) is available conjugated to agarose (sc-393194 AC), 500 µg/0.25 ml agarose in 1 ml, for IP; to HRP (sc-393194 HRP), 200 µg/ml, for WB, IHC(P) and ELISA; to either phycoerythrin (sc-393194 PE), fluorescein (sc-393194 FITC), Alexa Fluor® 488 (sc-393194 AF488), Alexa Fluor® 546 (sc-393194 AF546), Alexa Fluor® 594 (sc-393194 AF594) or Alexa Fluor® 647 (sc-393194 AF647), 200 µg/ml, for WB (RGB), IF, IHC(P) and FCM; and to either Alexa Fluor® 680 (sc-393194 AF680) or Alexa Fluor® 790 (sc-393194 AF790), 200 µg/ml, for Near-Infrared (NIR) WB, IF and FCM.

Alexa Fluor® is a trademark of Molecular Probes, Inc., Oregon, USA

APPLICATIONS

APIP (G-2) is recommended for detection of APIP of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for APIP siRNA (h): sc-61976, APIP siRNA (m): sc-61977, APIP shRNA Plasmid (h): sc-61976-SH, APIP shRNA Plasmid (m): sc-61977-SH, APIP shRNA (h) Lentiviral Particles: sc-61976-V and APIP shRNA (m) Lentiviral Particles: sc-61977-V.

Molecular Weight of APIP: 27 kDa.


Positive Controls: APIP (h): 293 Lysate: sc-111191, HEL 92.1.7 cell lysate: sc-2270 or HeLa whole cell lysate: sc-2200.

RECOMMENDED SUPPORT REAGENTS


To ensure optimal results, the following support reagents are recommended:

- 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.
- 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).
- 3) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

DATA



APIP (G-2): sc-393194. Western blot analysis of APIP expression in HeLa (A), HEL 92.1.7 (B) and F9 (C) whole cell lysates.



APIP (G-2): sc-393194. Western blot analysis of APIP expression in non-transfected: sc-110760 (A) and human APIP transfected: sc-111191 (B) 293 whole cell lysates.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.