

Cardiotrophin-1 siRNA (m): sc-39328

BACKGROUND

Cardiotrophin-1 (CT-1) is a member of the IL-6 family of cytokines, which signal through gp130 receptor complexes. gp130 complexes with several different receptor subunits to transmit signals from Cardiotrophin-1, IL-6, LIF, OSM, CNTF and IL-11. Cardiotrophin-1 binds to and activates the leukemia inhibitory factor (LIF) receptor/gp130 receptor complex and has been shown to induce hypertrophy in cardiac myocytes *in vitro*. Cardiotrophin-1, a secreted protein expressed at high levels in myocardium during cardiogenesis, has been shown to promote proliferation and survival of embryonic cardiomyocytes, suggesting a role for Cardiotrophin-1 in the activation of gp130 during cardiac development. Cardiotrophin-1 is highly expressed in heart, prostate, ovary and skeletal muscle. Lower levels of expression are seen in lung, kidney, pancreas, thymus, testis and small intestine.

REFERENCES

1. Hibi, M., et al. 1990. Molecular cloning and expression of an IL-6 signal transducer, gp130. *Cell* 63: 1149-1157.
2. Pennica, D., et al. 1995. Expression cloning of Cardiotrophin-1, a cytokine that induces cardiac myocyte hypertrophy. *Proc. Natl. Acad. Sci. USA* 92: 1142-1146.
3. Pennica, D., et al. 1996. Human Cardiotrophin-1: protein and gene structure, biological and binding activities, and chromosomal localization. *Cytokine* 8: 183-189.
4. Ishikawa, M., et al. 1996. cDNA cloning of rat Cardiotrophin-1 (CT-1): augmented expression of CT-1 gene in ventricle of genetically hypertensive rats. *Biochem. Biophys. Res. Commun.* 219: 377-381.
5. Wollert, K.C. and Chien, K.R. 1997. Cardiotrophin-1 and the role of gp130-dependent signaling pathways in cardiac growth and development. *J. Mol. Med.* 75: 492-501.
6. Klein, B. 1998. Update of gp130 cytokines in multiple myeloma. *Curr. Opin. Hematol.* 5: 186-191.

CHROMOSOMAL LOCATION

Genetic locus: Ctf1 (mouse) mapping to 7 F3.

PRODUCT

Cardiotrophin-1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Cardiotrophin-1 shRNA Plasmid (m): sc-39328-SH and Cardiotrophin-1 shRNA (m) Lentiviral Particles: sc-39328-V as alternate gene silencing products.

For independent verification of Cardiotrophin-1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-39328A, sc-39328B and sc-39328C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Cardiotrophin-1 siRNA (m) is recommended for the inhibition of Cardiotrophin-1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Cardiotrophin-1 (AN-B3): sc-9991 is recommended as a control antibody for monitoring of Cardiotrophin-1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Cardiotrophin-1 gene expression knockdown using RT-PCR Primer: Cardiotrophin-1 (m)-PR: sc-39328-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.