

IFN- β siRNA (h): sc-39603

BACKGROUND

The genes encoding type I interferons (IFNs), which include 14 IFN- α genes, 1 IFN- β gene, 1 IFN- ω (also known as IFN- α II1) gene, and a number of IFN- ω pseudogenes, are clustered on human chromosome 9p21.3. Interferons- α and - β are cytokines that are widely known to induce potent anti-viral activity. IFN- α and - β exert a variety of other biological effects, including anti-tumor and immunomodulatory activities and are increasingly used clinically to treat a range of malignancies, myelodysplasias and autoimmune diseases. IFN- ω is antigenically different from human IFN- α , IFN- β or IFN- γ , but is a component of natural mixtures of IFN species produced by virus-induced leukocytes or Burkitt's lymphoma cells. The Type I interferon receptor (IFN- α R) interacts with IFN- α , IFN- β and IFN- ω , and seems to be a multisubunit receptor.

REFERENCES

1. Adolf, G.R. 1987. Antigenic structure of human interferon ω 1 (interferon α II1): comparison with other human interferons. *J. Gen. Virol.* 68: 1669-1676.
2. Lim, J.K., et al. 1994. Intrinsic ligand binding properties of the human and bovine α -interferon receptors. *FEBS Lett.* 350: 281-286.
3. Hussain, M., et al. 1996. Identification of interferon- α 7, α 14, and α 21 variants in the genome of a large human population. *J. Interferon Cytokine Res.* 16: 853-859.
4. Mire-Sluis, A.R., et al. 1996. An anti-cytokine bioactivity assay for interferons- α , β and ω . *J. Immunol. Methods* 195: 55-61.
5. Cutrone, E.C., et al. 1997. Contributions of cloned type I interferon receptor subunits to differential ligand binding. *FEBS Lett.* 404: 197-202.
6. Vannucchi, S., et al. 2005. TRAIL is a key target in S-phase slowing-dependent apoptosis induced by interferon- β in cervical carcinoma cells. *Oncogene* 24: 2536-2546.
7. Siren, J., et al. 2005. IFN- α regulates TLR-dependent gene expression of IFN- α , IFN- β , IL-28 and IL-29. *J. Immunol.* 174: 1932-1937.
8. Molnarfi, N., et al. 2005. The production of IL-1 receptor antagonist in IFN- β -stimulated human monocytes depends on the activation of phosphatidylinositol 3-kinase but not of Stat1. *J. Immunol.* 174: 2974-2980.

CHROMOSOMAL LOCATION

Genetic locus: IFNB1 (human) mapping to 9p21.3.

PRODUCT

IFN- β siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see IFN- β shRNA Plasmid (h): sc-39603-SH and IFN- β shRNA (h) Lentiviral Particles: sc-39603-V as alternate gene silencing products.

For independent verification of IFN- β (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-39603A, sc-39603B and sc-39603C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

IFN- β siRNA (h) is recommended for the inhibition of IFN- β expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

IFN- β (A1): sc-53968 is recommended as a control antibody for monitoring of IFN- β gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor IFN- β gene expression knockdown using RT-PCR Primer: IFN- β (h)-PR: sc-39603-PR (20 μ l, 558 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.