

frizzled-2 siRNA (h): sc-39979

BACKGROUND

The frizzled gene, originally identified in *Drosophila melanogaster*, is involved in the development of tissue polarity. The mammalian homolog of frizzled as well as several secreted mammalian frizzled-related proteins (FRPs) have been described. The frizzled proteins contain seven transmembrane domains, a cysteine-rich domain in the extracellular region and a carboxy terminal Ser/Thr-xxx-Val motif. They function as receptors for Wnt and are generally coupled to G proteins. Expression of frizzled-2 can be observed in the fetal kidney and lung and in the adult ovary and colon. The Wnt/cGMP/Ca²⁺ pathway is mediated by frizzled-2. It binds Wnt proteins and signals by activating the release of stored calcium. Expression of frizzled-2 is regulated by Angiotensin II. Activated frizzled-2 suppresses the activity of protein kinase G, and activates NFAT-dependent transcription, the phosphatidylinositol pathway and calcium sensitive enzymes.

REFERENCES

- Zhao, Z., et al. 1995. A human homologue of the *Drosophila* polarity gene frizzled has been identified and mapped to 17q21.1. *Genomics* 27: 370-373.
- Sagara, N., et al. 1998. Molecular cloning, differential expression, and chromosomal localization of human frizzled-1, frizzled-2, and frizzled-7. *Biochem. Biophys. Res. Commun.* 252: 117-122.
- Malbon, C.C., et al. 2004. Frizzleds: new members of the superfamily of G protein-coupled receptors. *Front. Biosci.* 9: 1048-1058.
- Castoldi, G., et al. 2005. Angiotensin II modulates frizzled-2 receptor expression in rat vascular smooth muscle cells. *Clin. Sci.* 108: 523-530.
- Rodriguez, J., et al. 2005. SFRP1 regulates the growth of retinal ganglion cell axons through the Fz2 receptor. *Nat. Neurosci.* 8: 1301-1309.
- Jiang, F., et al. 2006. Gene expression profile of quiescent and activated rat hepatic stellate cells implicates Wnt signaling pathway in activation. *J. Hepatol.* 45: 401-409.
- Ma, L., et al. 2006. Suppression of cyclic GMP-dependent protein kinase is essential to the Wnt/cGMP/Ca²⁺ pathway. *J. Biol. Chem.* 281: 30990-31001.
- Harris, K.E., et al. 2007. Different Wnt signals act through the frizzled and RYK receptors during *Drosophila* salivary gland migration. *Development* 134: 2017-2025.
- Hurson, C.J., et al. 2007. Gene expression analysis in human osteoblasts exposed to dexamethasone identifies altered developmental pathways as putative drivers of osteoporosis. *BMC Musculoskelet. Disord.* 8: 12.

CHROMOSOMAL LOCATION

Genetic locus: FZD2 (human) mapping to 17q21.31.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

frizzled-2 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see frizzled-2 shRNA Plasmid (h): sc-39979-SH and frizzled-2 shRNA (h) Lentiviral Particles: sc-39979-V as alternate gene silencing products.

For independent verification of frizzled-2 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-39979A, sc-39979B and sc-39979C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

frizzled-2 siRNA (h) is recommended for the inhibition of frizzled-2 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor frizzled-2 gene expression knockdown using RT-PCR Primer: frizzled-2 (h)-PR: sc-39979-PR (20 μ l, 554 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.