

# CRALBP siRNA (h): sc-40428

## BACKGROUND

11-*cis*-retinal, the universal chromophore of the vertebrate retina, is coupled to opsins in both rod and cone photoreceptor cells and is photoisomerized to all-*trans*-retinal by light. This conversion is inhibited when 11-*cis*-retinol is in a complex with cellular retinaldehyde-binding protein (CRALBP). CRALBP may play a role in the vertebrate visual process as a substrate-routing protein, influencing the enzymatic partitioning of 11-*cis*-retinol at a key branch point in the visual cycle. Human CRALBP maps to chromosome 15q26.1 and encodes a 316 amino acid protein. CRALBP is not expressed in photoreceptors and is abundant in the retinal pigment epithelium (RPE) and Muller cells of the neuroretina, where it carries 11-*cis*-retinol and 11-*cis*-retinaldehyde. Mutations in the human CRALBP gene cause retinal pathology and delayed dark adaptation. CRALBP knockout mice have a delayed response in rhodopsin regeneration, 11-*cis*-retinal production and dark adaptation after illumination.

## REFERENCES

- Crabb, J.W., et al. 1988. Cloning of the cDNAs encoding the cellular retinaldehyde-binding protein from bovine and human retina and comparison of the protein structures. *J. Biol. Chem.* 263: 18688-18692.
- Sparkes, R.S., et al. 1992. Assignment of the gene (RLBP1) for cellular retinaldehyde-binding protein (CRALBP) to human chromosome 15q26 and mouse chromosome 7. *Genomics* 12: 58-62.
- Intres, R., et al. 1994. Molecular cloning and structural analysis of the human gene encoding cellular retinaldehyde-binding protein. *J. Biol. Chem.* 269: 25411-25418.
- McBee, J.K., et al. 2001. Isomerization of 11-*cis*-retinoids to all-*trans*-retinoids *in vitro* and *in vivo*. *J. Biol. Chem.* 276: 48483-48493.
- McBee, J.K., et al. 2001. Visual cycle impairment in cellular retinaldehyde binding protein (CRALBP) knockout mice results in delayed dark adaptation. *Neuron* 29: 739-748.
- Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2001. Johns Hopkins University, Baltimore, MD. MIM Number: 180090. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

## CHROMOSOMAL LOCATION

Genetic locus: RLBP1 (human) mapping to 15q26.1.

## PRODUCT

CRALBP siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CRALBP shRNA Plasmid (h): sc-40428-SH and CRALBP shRNA (h) Lentiviral Particles: sc-40428-V as alternate gene silencing products.

For independent verification of CRALBP (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-40428A, sc-40428B and sc-40428C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

CRALBP siRNA (h) is recommended for the inhibition of CRALBP expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

CRALBP (G-9): sc-376082 is recommended as a control antibody for monitoring of CRALBP gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CRALBP gene expression knockdown using RT-PCR Primer: CRALBP (h)-PR: sc-40428-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.