

β-1,4-Gal-T4 siRNA (m): sc-40621

BACKGROUND

Enzymatic glycosylation of proteins and lipids is an important biological process. A large number of glycosyltransferases synthesize a wide variety of glycoconjugates. A novel putative member of the human UDP-galactose: β-N-acetylglucosamine β-1,4-galactosyltransferase family, designated β-1,4-Gal-T4, encodes a type II membrane protein which has significant sequence similarity to other β-1,4-galactosyltransferases. β-1,4-Gal-T4 catalyzes glycosylation of glycolipids with terminal β-GlcNAc. Unlike β-1,4-Gal-T1, -T2 and -T3, β-1,4-Gal-T4 does not transfer galactose to asialo-agalacto-fetuin, asialo-agacto-transferrin or ovalbumin. β-1,4-Gal-T4 has a very restricted pattern of tissue expression. β-1,4-Gal-T4 is localized to two subcellular compartments, the Golgi complex, where it participates in cellular glycosylation, and the plasma membrane, where it functions as a receptor for oligosaccharide ligands on opposing cells or in the extracellular matrix.

REFERENCES

1. Shur, B.D. 1984. The receptor function of galactosyltransferase during cellular interactions. *Mol. Cell. Biochem.* 61: 143-158.
2. Shur, B.D. 1986. The receptor function of galactosyltransferase during mammalian fertilization. *Adv. Exp. Med. Biol.* 207: 79-93.
3. Strous, G.J. 1986. Golgi and secreted galactosyltransferase. *CRC Crit. Rev. Biochem.* 21: 119-151.
4. Amado, M., et al. 1998. A family of human β3-galactosyltransferases. Characterization of four members of a UDP-galactose:β-N-acetylglucosamine/β-nacetyl-galactosamine β-1,3-galactosyltransferase family. *J. Biol. Chem.* 273: 12770-12778.
5. Amado, M., et al. 1999. Identification and characterization of large galactosyltransferase gene families: galactosyltransferases for all functions. *Biochim. Biophys. Acta* 1473: 35-53.
6. Online Mendelian Inheritance in Man, OMIM™. 2001. Johns Hopkins University, Baltimore, MD. MIM Number: 603093. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
7. LocusLink Report (LocusID: 2683). <http://www.ncbi.nlm.nih.gov/LocusLink/>

CHROMOSOMAL LOCATION

Genetic locus: B4galt4 (mouse) mapping to 16 B4.

PRODUCT

β-1,4-Gal-T4 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μM solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see β-1,4-Gal-T4 shRNA Plasmid (m): sc-40621-SH and β-1,4-Gal-T4 shRNA (m) Lentiviral Particles: sc-40621-V as alternate gene silencing products.

For independent verification of β-1,4-Gal-T4 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-40621A, sc-40621B and sc-40621C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μl of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μl of RNase-free water makes a 10 μM solution in a 10 μM Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

β-1,4-Gal-T4 siRNA (m) is recommended for the inhibition of β-1,4-Gal-T4 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μM in 66 μl. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor β-1,4-Gal-T4 gene expression knockdown using RT-PCR Primer: β-1,4-Gal-T4 (m)-PR: sc-40621-PR (20 μl, 597 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.