



# Selenoprotein W siRNA (h): sc-40932

## BACKGROUND

Selenium is an essential trace element that is incorporated as selenocysteine into the primary structure of selenoproteins. Nutritional deficiency of selenium decreases selenoprotein concentrations and leads to pathologic conditions. Most of the known selenoproteins are members of the Glutathione peroxidase or iodothyronine deiodinase families. The Selenoprotein N glycoprotein localizes to the endoplasmic reticulum (ER) and contains selenocysteine at its active site. There are two isoforms associated with Selenoprotein N: isoform 1, the full-length transcript; and isoform 2, which lacks exon 3. Selenoprotein N is primarily expressed in skeletal muscle, brain, lung and placenta, but isoform 2 can also be detected in heart and stomach tissues. Mutations in SEP1, the gene encoding for Selenoprotein, cause multimicore disease and rigid spine muscular dystrophy.

## REFERENCES

1. Keshan Disease Research Group. Observations on effect of sodium selenite in prevention of Keshan disease. 1979. Chin. Med. J. 92: 471-476.
2. Vendeland, S.C., Beilstein, M.A., Chen, C.L., Jensen, O.N., Barofsky, E. and Whanger, P.D. 1993. Purification and properties of Selenoprotein W from rat muscle. J. Biol. Chem. 268: 17103-17107.
3. Vendeland, S.C., Beilstein, M.A., Yeh, J.Y., Ream, W. and Whanger, P.D. 1995. Rat skeletal muscle Selenoprotein W: cDNA clone and mRNA modulation by dietary selenium. Proc. Natl. Acad. Sci. USA 92: 8749-8753.
4. Hill, K.E., Dasouki, M., Phillips, J.A., III. and Burk, R.F. 1996. Human Selenoprotein P gene maps to 5q31. Genomics 36: 550-551.
5. Gu, Q.P., Beilstein, M.A., Vendeland, S.C., Lugade, A., Ream, W. and Whanger, P.D. 1997. Conserved features of selenocysteine insertion sequence (SECIS) elements in Selenoprotein W cDNAs from five species. Gene 193: 187-196.
6. Gu, Q.P., Sun, Y., Ream, L.W. and Whanger, P.D. 2000. Selenoprotein W accumulates primarily in primate skeletal muscle, heart, brain and tongue. Mol. Cell. Biochem. 204: 49-56.
7. Jeong, D., Kim, T.S., Chung, Y.W., Lee, B.J. and Kim, I.Y. 2002. Selenoprotein W is a glutathione-dependent antioxidant *in vivo*. FEBS Lett. 517: 225-228.
8. Jeong, D.W., Kim, E.H., Kim, T.S., Chung, Y.W., Kim, H. and Kim, I.Y. 2004. Different distributions of selenoprotein W and thioredoxin during postnatal brain development and embryogenesis. Mol. Cells 17: 156-159.
9. Loflin, J., Lopez, N., Whanger, P.D. and Kioussi, C. 2006. Selenoprotein W during development and oxidative stress. J. Inorg. Biochem. 100: 1679-1684.

## CHROMOSOMAL LOCATION

Genetic locus: SEP1 (human) mapping to 19q13.33.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

Selenoprotein W siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Selenoprotein W shRNA Plasmid (h): sc-40932-SH and Selenoprotein W shRNA (h) Lentiviral Particles: sc-40932-V as alternate gene silencing products.

For independent verification of Selenoprotein W (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-40932A, sc-40932B and sc-40932C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

Selenoprotein W siRNA (h) is recommended for the inhibition of Selenoprotein W expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Selenoprotein W gene expression knockdown using RT-PCR Primer: Selenoprotein W (h)-PR: sc-40932-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## SELECT PRODUCT CITATIONS

1. Ren, B., Situ, J., Huang, X., Tan, Q., Xiao, S., Li, N., Tian, J., Du, X., Ni, J. and Liu, Q. 2024. Selenoprotein W modulates Tau homeostasis in an Alzheimer's disease mouse model. Commun. Biol. 7: 872.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.