

# ABCD2 siRNA (h): sc-41145

## BACKGROUND

The peroxisomal membrane contains several ATP-binding cassette (ABC) transporters, designated ABCD1-4, that are known to be present in the human peroxisome membrane. All four proteins are ABC half-transporters, which dimerize to form an active transporter. A mutation in the ABCD1 gene causes X-linked adrenoleukodystrophy (X-ALD), a peroxisomal disorder which affects lipid storage. ABCD2 in mouse is expressed at high levels in the brain and adrenal organs, which are adversely affected in X-ALD. The peroxisomal membrane comprises two quantitatively major proteins, PMP22 and ABCD3. ABCD3 is associated with irregularly shaped vesicles which may be defective peroxisomes or peroxisome precursors. ABCD1 localizes to peroxisomes. ABCB7 is a half-transporter involved in the transport of heme from the mitochondria to the cytosol.

## REFERENCES

1. Gartner, J., et al. 1992. Characterization and localization of the human 70 kDa peroxisomal membrane protein (PMP70) gene. *Am. J. Hum. Genet.* 51: 168.
2. Lombard-Platet, G., et al. 1996. A close relative of the adrenoleukodystrophy (ALD) gene codes for a peroxisomal protein with a specific expression pattern. *Proc. Natl. Acad. Sci. USA* 93: 1265-1269.
3. Shani, N., et al. 1997. Identification of a fourth half ABC transporter in the human peroxisomal membrane. *Hum. Mol. Genet.* 6: 1925-1931.
4. Moser, H.W. 1997. Adrenoleukodystrophy: phenotype, genetics, pathogenesis and therapy. *Brain* 120: 1485-1508.
5. Savary, S., et al. 1997. Chromosomal localization of the adrenoleukodystrophy-related gene in man and mice. *Eur. J. Hum. Genet.* 5: 99-101.
6. Holzinger, A., et al. 1998. Genomic organization and chromosomal localization of the human peroxisomal membrane protein-1-like protein (PMP1-L) gene encoding a peroxisomal ABC transporter. *FEBS Lett.* 426: 238-242.
7. Shimada, Y., et al. 1998. Cloning and chromosomal mapping of a novel ABC transporter gene (hABC7), a candidate for X-linked sideroblastic anemia with spinocerebellar ataxia. *J. Hum. Genet.* 43: 115-122.
8. LocusLink Report (LocusID: 300100). <http://www.ncbi.nlm.nih.gov/LocusLink/>

## CHROMOSOMAL CHROMOSOMAL LOCATION

Genetic locus: ABCD2 (human) mapping to 12q12.

## PRODUCT

ABCD2 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ABCD2 shRNA Plasmid (h): sc-41145-SH and ABCD2 shRNA (h) Lentiviral Particles: sc-41145-V as alternate gene silencing products.

For independent verification of ABCD2 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-41145A, sc-41145B and sc-41145C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

ABCD2 siRNA (h) is recommended for the inhibition of ABCD2 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ABCD2 gene expression knockdown using RT-PCR Primer: ABCD2 (h)-PR: sc-41145-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## SELECT PRODUCT CITATIONS

1. Singh, J., et al. 2011. HDAC inhibitor SAHA normalizes the levels of VLCFAs in human skin fibroblasts from X-ALD patients and downregulates the expression of proinflammatory cytokines in ABCD1/2-silenced mouse astrocytes. *J. Lipid Res.* 52: 2056-2069.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.