

ILBP siRNA (h): sc-41241

BACKGROUND

The fatty acid binding protein (FABP) family of cytoplasmic hydrophobic ligand binding proteins influence lipid metabolism by binding and transporting long-chain fatty acids. Ileal lipid binding protein (ILBP) is a cytosolic ileocyte FABP that binds to both bile acids and fatty acids thereby mediating active uptake of bile acid in the ileum. Transport of bile acids from the liver is essential for the solubilization and transport of dietary lipids. ILBP contains ten antiparallel β strands arranged in two nearly orthogonal β sheets (β clam shell), covered on one side by two short, nearly parallel α helices. Binding of fatty acids or bile acids to ILBP alters the side-chain proton resonances of amino acids within the protein cavity and increases the affinity of ILBP for bile acids; bile acid binding to ILBP is a positive-feedback regulation mechanism. The human ILBP gene maps to position 5q33.3, with transcript being abundant in the small intestine.

REFERENCES

- Oelkers, P. and Dawson, P.A. 1995. Cloning and chromosomal localization of the human ileal lipid-binding protein. *Biochim. Biophys. Acta* 1257: 199-202.
- Lucke, C., et al. 1996. Flexibility is a likely determinant of binding specificity in the case of ileal lipid binding protein. *Structure* 4: 785-800.
- Borchers, T., et al. 1997. Heart-type fatty acid binding protein-involvement in growth inhibition and differentiation. *Prostaglandins Leukot. Essent. Fatty Acids* 57: 77-84.
- Kramer, W., et al. 1998. Intestinal absorption of bile acids: paradoxical behaviour of the 14 kDa ileal lipid-binding protein in differential photo-affinity labelling. *Biochem. J.* 333: 335-341.
- Lucke, C., et al. 1999. A comparative study of the backbone dynamics of two closely related lipid binding proteins: bovine heart fatty acid binding protein and porcine ileal lipid binding protein. *Mol. Cell. Biochem.* 192: 109-121.
- Grober, J., et al. 1999. Identification of a bile acid-responsive element in the human ileal bile acid-binding protein gene. Involvement of the farnesoid X receptor/9-*cis*-retinoic acid receptor heterodimer. *J. Biol. Chem.* 274: 29749-29754.

CHROMOSOMAL LOCATION

Genetic locus: FABP6 (human) mapping to 5q33.3.

PRODUCT

ILBP siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ILBP shRNA Plasmid (h): sc-41241-SH and ILBP shRNA (h) Lentiviral Particles: sc-41241-V as alternate gene silencing products.

For independent verification of ILBP (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-41241A, sc-41241B and sc-41241C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ILBP siRNA (h) is recommended for the inhibition of ILBP expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

ILBP (E-9): sc-515609 is recommended as a control antibody for monitoring of ILBP gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz MarkerTM Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ILBP gene expression knockdown using RT-PCR Primer: ILBP (h)-PR: sc-41241-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.