

# ADH $\gamma$ siRNA (h) : sc-41440

## BACKGROUND

The alcohol dehydrogenase family of proteins metabolize a wide variety of substrates, including ethanol, retinol, other aliphatic alcohols, hydroxysteroids, and lipid peroxidation products. Class I alcohol dehydrogenase, consisting of several homo- and heterodimers of  $\alpha$ ,  $\beta$ , and  $\gamma$  subunits, exhibits high activity for ethanol oxidation and plays a major role in ethanol catabolism. Three genes encoding  $\alpha$  (ADH1A),  $\beta$  (ADH1B) and  $\gamma$  (ADH1C) subunits are tandemly organized on chromosome 4q23 as a gene cluster. The  $\alpha$  form of ADH is monomeric and predominant in fetal and infant livers, becoming less active in gestation and only weakly active during adulthood. The genes encoding  $\beta$  and  $\gamma$  subunits, however, are polymorphic and strongly expressed in adult livers. With the coenzyme NAD, ADH catalyzes the reversible conversion of organic alcohols to ketones or aldehydes. The physiologic function for ADH in the liver is the removal of ethanol formed by microorganisms in the intestinal tract.

## REFERENCES

1. Smith, M., et al. 1973. Studies on the subunit structure and molecular size of the human dehydrogenase isozymes determined by the different loci, ADH1, ADH2, and ADH3. *Ann. Hum. Genet.* 36: 401-414.
2. Smith, M., et al. 1984. Derivation of probes for molecular genetic analysis of human class I alcohol dehydrogenase (ADH), a polymorphic gene family on chromosome 4. *Am. J. Hum. Genet.* 36: 153S.
3. Tsukahara, M. and Yoshida, A. 1989. Chromosomal assignment of the alcohol dehydrogenase cluster locus to human chromosome 4q21-23 by *in situ* hybridization. *Genomics* 4: 218-220.
4. Yasunami, M., et al. 1989. The organization of human class I alcohol dehydrogenase gene cluster. *Cytogenet. Cell Genet.* 51: 1113.
5. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 103700. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Jelski, W., et al. 2007. Alcohol dehydrogenase (ADH) isoenzyme activity in the sera of patients with *Helicobacter pylori* infection. *Dig. Dis. Sci.* 52: 1513-1516.

## CHROMOSOMAL LOCATION

Genetic locus: ADH1C (human) mapping to 4q23.

## PRODUCT

ADH $\gamma$  siRNA (h) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ADH $\gamma$  shRNA Plasmid (h): sc-41440-SH and ADH $\gamma$  shRNA (h) Lentiviral Particles: sc-41440-V as alternate gene silencing products.

For independent verification of ADH $\gamma$  (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-41440A and sc-41440B.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

ADH $\gamma$  siRNA (h) is recommended for the inhibition of ADH $\gamma$  expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

ADH (G-7): sc-133207 is recommended as a control antibody for monitoring of ADH $\gamma$  gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ADH $\gamma$  gene expression knockdown using RT-PCR Primer: ADH $\gamma$  (h)-PR: sc-41440-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.