

GlyT1 siRNA (h): sc-41974

BACKGROUND

Na⁺/Cl⁻ dependent neurotransmitter transporters are a superfamily of transmembrane proteins that contain 12 membrane spanning regions. Specifically, the highly hydrophobic Na⁺/Cl⁻ dependent glycine transporters (GlyT) are crucial for the termination of neurotransmission at glycinergic synapses. Two different GlyT genes encode GlyT2 and GlyT1, which exists as two isoforms produced by alternative splicing of the same gene located on human chromosome 1p34.1. The GlyT1 gene may be an early marker of neural development and encodes glia-specific transporter proteins. Although GlyT1 and GlyT2 are both expressed in the brain and spinal cord, each shows a unique pattern of expression. GlyT1 is found only in the white matter of the CNS, whereas GlyT2 is found in the gray matter of the CNS as well as in macrophages and mast cells in peripheral tissues. The anatomic distribution of GlyT2 mRNA suggests that glycine may act as a supraspinal neurotransmitter and may function as a chemical messenger outside the CNS.

REFERENCES

1. Liu, Q.R., Nelson, H., Mandiyan, S., Lopez-Corcuera, B. and Nelson, N. 1992. Cloning and expression of a glycine transporter from mouse brain. *FEBS Lett.* 305: 110-114.
2. Borowsky, B., Mezey, E. and Hoffman, B.J. 1993. Two glycine transporter variants with distinct localization in the CNS and peripheral tissues are encoded by a common gene. *Neuron* 10: 851-863.
3. Kim, K.M., Kingsmore, S.F., Han, H., Yang-Feng, T.L., Godinot, N., Seldin, M.F., Caron, M.G. and Giros, B. 1994. Cloning of the human glycine transporter type 1: molecular and pharmacological characterization of novel isoform variants and chromosomal localization of the gene in the human and mouse genomes. *Mol. Pharmacol.* 45: 608-617.
4. Adams, R.H., Sato, K., Shimada, S., Tohyama, M., Puschel, A.W. and Betz, H. 1995. Gene structure and glial expression of the glycine transporter GlyT1 in embryonic and adult rodents. *J. Neurosci.* 15: 2524-2532.
5. Evans, J., Herdon, H., Cairns, W., O'Brien, E., Chapman, C., Terrett, J. and Gloger, I. 1999. Cloning, functional characterisation and population analysis of a variant form of the human glycine type 2 transporter. *FEBS Lett.* 463: 301-306.
6. Horiuchi, M., Nicke, A., Gomez, J., Ashcraft, A., Schmalzing, G. and Betz, H. 2001. Surface-localized glycine transporters 1 and 2 function as monomeric proteins in *Xenopus* oocytes. *Proc. Natl. Acad. Sci. USA* 98: 1448-1453.

CHROMOSOMAL LOCATION

Genetic locus: SLC6A9 (human) mapping to 1p34.1.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

GlyT1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GlyT1 shRNA Plasmid (h): sc-41974-SH and GlyT1 shRNA (h) Lentiviral Particles: sc-41974-V as alternate gene silencing products.

For independent verification of GlyT1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-41974A, sc-41974B and sc-41974C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GlyT1 siRNA (h) is recommended for the inhibition of GlyT1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GlyT1 gene expression knockdown using RT-PCR Primer: GlyT1 (h)-PR: sc-41974-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.