

TRPC6 siRNA (h): sc-42672

BACKGROUND

Transient receptor potential cation (TRPC) channels are a superfamily of six transmembrane segment-spanning, gated cation channels. TRPC subtypes mediate store-operated Ca^{2+} entry, a process involving Ca^{2+} influx and replenishment of Ca^{2+} stores formerly emptied through the action of inositol 1,4,5-trisphosphate production and other Ca^{2+} mobilizing agents. TRPC ion channels influence calcium-depletion induced calcium influx processes in response to chemo-, mechano- and osmoregulatory events. Human TRPC6 protein is a 931 amino acid cation channel that is predominantly expressed in placenta, spleen, lung, small intestine and ovary. Activated by diacylglycerol (DAG), TRPC6 comprises the α_1 -adrenoceptor-activated Ca^{2+} -permeable cation channel. The gene encoding human TRPC6 maps to chromosome 11q22.1.

CHROMOSOMAL LOCATION

Genetic locus: TRPC6 (human) mapping to 11q22.1.

PRODUCT

TRPC6 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μM solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TRPC6 shRNA Plasmid (h): sc-42672-SH and TRPC6 shRNA (h) Lentiviral Particles: sc-42672-V as alternate gene silencing products.

For independent verification of TRPC6 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-42672A, sc-42672B and sc-42672C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20°C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20°C , avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μl of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μl of RNase-free water makes a 10 μM solution in a 10 μM Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

TRPC6 siRNA (h) is recommended for the inhibition of TRPC6 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μM in 66 μl . Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

TRPC6 (B-10): sc-515837 is recommended as a control antibody for monitoring of TRPC6 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TRPC6 gene expression knockdown using RT-PCR Primer: TRPC6 (h)-PR: sc-42672-PR (20 μl , 516 bp). Annealing temperature for the primers should be $55-60^\circ\text{C}$ and the extension temperature should be $68-72^\circ\text{C}$.

SELECT PRODUCT CITATIONS

- Antigny, F., et al. 2011. Transient receptor potential canonical channel 6 links Ca^{2+} mishandling to cystic fibrosis transmembrane conductance regulator channel dysfunction in cystic fibrosis. *Am. J. Respir. Cell Mol. Biol.* 44: 83-90.
- Liao, C., et al. 2012. The upregulation of TRPC6 contributes to Ca^{2+} signaling and Actin assembly in human mesangial cells after chronic hypoxia. *Biochem. Biophys. Res. Commun.* 421: 750-756.
- Meng, K., et al. 2014. Calcium sensing receptor modulates extracellular calcium entry and proliferation via TRPC3/6 channels in cultured human mesangial cells. *PLoS ONE* 9: e98777.
- Zhao, B., et al. 2015. The role of TRPC6 in oxidative stress-induced podocyte ischemic injury. *Biochem. Biophys. Res. Commun.* 461: 413-420.
- Wen, L., et al. 2016. Regulation of multi-drug resistance in hepatocellular carcinoma cells is TRPC6/calcium dependent. *Sci. Rep.* 6: 23269.
- Zhang, L., et al. 2017. Calcium-sensing receptor stimulation in cultured glomerular podocytes induces TRPC6-dependent calcium entry and Rho A activation. *Cell. Physiol. Biochem.* 43: 1777-1789.
- Wang, Q., et al. 2019. Group I metabotropic glutamate receptor activation induces TRPC6-dependent calcium influx and Rho A activation in cultured human kidney podocytes. *Biochem. Biophys. Res. Commun.* 511: 374-380.
- Wang, D., et al. 2020. The role of TRPC6 in α_1 -AR activation-induced calcium signal changes in human podocytes. *Ann. Palliat. Med.* 9: 1596-1605.
- Hsu, W.L., et al. 2020. Nociceptive transient receptor potential canonical 7 (TRPC7) mediates aging-associated tumorigenesis induced by ultraviolet B. *Aging Cell* 19: e13075.
- Yart, L., et al. 2022. Dual effect of nifedipine on pregnant human myometrium contractility: implication of TRPC1. *J. Cell. Physiol.* 237: 1980-1991.
- Yu, H., et al. 2022. TRPC6 mediates high glucose-induced mitochondrial fission through activation of CDK5 in cultured human podocytes. *Front. Physiol.* 13: 984760.

RESEARCH USE

For research use only, not for use in diagnostic procedures.