

Myosin VIIa siRNA (m): sc-43224

BACKGROUND

Myosins are molecular motors that move along filamentous Actin and influence cellular movements such as phagocytosis. There are seven classes of myosins in vertebrates, including Myosin II, and six unconventional Myosin classes, designated I, V, VI, VII, IX and X. Myosin VIIa is a plus end-directed motor that influences cilia formation and cell adhesion. Mutations in the human Myosin VIIa gene correlate with Usher syndrome, a disease characterized by congenital sensorineural deafness, vestibular dysfunction and retinitis pigmentosa.

REFERENCES

1. Weil, D., et al. 1995. Defective Myosin VIIa gene responsible for Usher syndrome type 1B. *Nature* 374: 60-61.
2. Weil, D., et al. 1996. Human Myosin VIIA responsible for the Usher 1B syndrome: a predicted membrane-associated motor protein expressed in developing sensory epithelia. *Proc. Natl. Acad. Sci. USA* 93: 3232-3237.
3. Weil, D., et al. 1997. The autosomal recessive isolated deafness, DFNB2, and the Usher 1B syndrome are allelic defects of the Myosin VIIa gene. *Nat. Genet.* 16: 191-193.
4. Maniak, M. 2001. Cell adhesion: ushering in a new understanding of Myosin VII. *Curr. Biol.* 11: R315-R317.
5. Tuxworth, R.I., et al. 2001. A role for Myosin VII in dynamic cell adhesion. *Curr. Biol.* 11: 318-329.
6. Online Mendelian Inheritance in Man, OMIM[™]. 2003. Johns Hopkins University, Baltimore, MD. MIM Number: 276903. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
7. LocusLink Report (LocusID: 4647). <http://www.ncbi.nlm.nih.gov/LocusLink/>

CHROMOSOMAL LOCATION

Genetic locus: Myo7a (mouse) mapping to 7 E2.

PRODUCT

Myosin VIIa siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Myosin VIIa shRNA Plasmid (m): sc-43224-SH and Myosin VIIa shRNA (m) Lentiviral Particles: sc-43224-V as alternate gene silencing products.

For independent verification of Myosin VIIa (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-43224A, sc-43224B and sc-43224C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Myosin VIIa siRNA (m) is recommended for the inhibition of Myosin VIIa expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Myosin VIIa (C-5): sc-74516 is recommended as a control antibody for monitoring of Myosin VIIa gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Myosin VIIa gene expression knockdown using RT-PCR Primer: Myosin VIIa (m)-PR: sc-43224-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.