

dysferlin siRNA (h): sc-43739

BACKGROUND

Dysferlin is a muscle-specific protein that is essential for normal muscle function and development. Mutations in the human dysferlin gene, DYSF, which maps to chromosome 2p13.2, are associated with limb girdle muscular dystrophy-2B (LGMD-2B) and a related, adult-onset, distal dystrophy known as Miyoshi myopathy (MM). Dysferlin localizes to the muscle fiber membrane, but is absent in MM and LGMD-2B muscle. Dysferlin is detected in 5-6 week embryos, when limbs begin to form regional differentiation. Although it is not essential for initial myogenesis, dysferlin appears to be critical for sustained normal function in mature muscle. It has been suggested that the absence of dysferlin during development gives rise to the disease phenotype in adulthood. Identical mutations in the dysferlin gene can produce more than one myopathy phenotype, indicating that additional genes and/or other factors are also involved in the clinical phenotype. The DYSF gene has no homology to any other known mammalian gene, but the protein product is related to the spermatogenesis factor fer-1 of *Caenorhabditis elegans*. The name "dysferlin" combines the role of the gene in producing muscular dystrophy with its homology to *C. elegans*.

REFERENCES

1. Bashir, R., et al. 1994. A gene for autosomal recessive limb-girdle muscular dystrophy maps to chromosome 2p. *Hum. Mol. Genet.* 3: 455-457.
2. Liu, J., et al. 1998. Dysferlin, a novel skeletal muscle gene, is mutated in Miyoshi myopathy and limb girdle muscular dystrophy. *Nat. Genet.* 20: 31-36.
3. Matsuda, C., et al. 1999. Dysferlin is a surface membrane-associated protein that is absent in Miyoshi myopathy. *Neurology* 53: 1119-1122.
4. Anderson, L.V., et al. 1999. Dysferlin is a plasma membrane protein and is expressed early in human development. *Hum. Mol. Genet.* 8: 855-861.
5. Weiler, T., et al. 1999. Identical mutation in patients with limb girdle muscular dystrophy type 2B or Miyoshi myopathy suggests a role for modifier gene(s). *Hum. Mol. Genet.* 8: 871-877.
6. LocusLink Report (LocusID: 829). <http://www.ncbi.nlm.nih.gov/LocusLink>

CHROMOSOMAL LOCATION

Genetic locus: DYSF (human) mapping to 2p13.2.

PRODUCT

dysferlin siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see dysferlin shRNA Plasmid (h): sc-43739-SH and dysferlin shRNA (h) Lentiviral Particles: sc-43739-V as alternate gene silencing products.

For independent verification of dysferlin (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-43739A, sc-43739B and sc-43739C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

dysferlin siRNA (h) is recommended for the inhibition of dysferlin expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

dysferlin (C-11): sc-398905 is recommended as a control antibody for monitoring of dysferlin gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor dysferlin gene expression knockdown using RT-PCR Primer: dysferlin (h)-PR: sc-43739-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

1. Azar, C., et al. 2021. RNA-Seq identifies genes whose proteins are upregulated during syncytia development in murine C2C12 myoblasts and human BeWo trophoblasts. *Physiol. Rep.* 9: e14671.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.