



LIF (hBA-FL): sc-4377

BACKGROUND

Embryonic stem (ES) cells are the focus of much research and represent great therapeutic potential as they can be propagated indefinitely in an undifferentiated state while possessing the ability to differentiate into all embryonic germ layers (endoderm, ectoderm and mesoderm) both *in vivo* and *in vitro*. LIF (leukemia inhibitory factor), also known as MLPLI (melanoma-derived LPL inhibitor), HILDA, DIA or CDF, is a 202 amino acid secreted protein and lymphoid factor that participates in the maintenance of ES cell pluripotency by suppressing spontaneous ES cell differentiation. Secreted LIF precursor is further processed into a biologically active glycoprotein. Expressed by a wide variety of cells including activated T lymphocytes, monocytes, mast cells and neuronal cells, LIF is suggested to promote survival and growth of axons *in vitro* and is involved in immune tolerance at the maternal-fetal interface. LIF may also participate in fat and bone metabolism and regulate epithelial conversion during kidney development.

REFERENCES

1. Gough, N.M., et al. 1988. Molecular cloning and expression of the human homologue of the murine gene encoding myeloid leukemia-inhibitory factor. *Proc. Natl. Acad. Sci. USA* 85: 2623-2627.
2. Patterson, P.H. 1994. Leukemia inhibitory factor, a cytokine at the interface between neurobiology and immunology. *Proc. Natl. Acad. Sci. USA* 91: 7833-7835.
3. Barasch, J., et al. 1999. Mesenchymal to epithelial conversion in rat metanephros is induced by LIF. *Cell* 99: 377-386.
4. Pera, M.F., et al. 2000. Human embryonic stem cells. *J. Cell Sci.* 113: 5-10.
5. Hu, W., et al. 2007. p53 regulates maternal reproduction through LIF. *Nature* 450: 721-724.
6. Niwa, H., et al. 2009. A parallel circuit of LIF signalling pathways maintains pluripotency of mouse ES cells. *Nature* 460: 118-122.
7. Online Mendelian Inheritance in Man, OMIM™. 2009. Johns Hopkins University, Baltimore, MD. MIM Number: 159540. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
8. Xu, J., et al. 2010. Role of leukaemia inhibitory factor in the induction of pluripotent stem cells in mice. *Cell Biol. Int.* 34: 791-797.
9. Fukunaga, N., et al. 2010. Leukemia inhibitory factor (LIF) enhances germ cell differentiation from primate embryonic stem cells. *Cell. Reprogram.* 12: 369-376.

SOURCE

LIF (hBA-FL) is an antibody of leukemia inhibitory factor (LIF) of human origin.

PRODUCT

LIF (hBA-FL) is purified from bacterial lysates (> 98%); supplied as 50 µg.

STORAGE

Store desiccated at -20° C; stable for one year from the date of shipment.

APPLICATIONS

LIF (hBA-FL) is a potent inhibitor of spontaneous embryonic mouse and human stem cell differentiation.

Molecular Weight of LIF precursor: 22 kDa.

Molecular Weight of mature glycosylated LIF: 40-45 kDa.

RECONSTITUTION

In order to avoid freeze/thaw damaging of the active protein, dilute protein when first used to desired working concentration. Either a sterile filtered standard buffer (such as 50mM TRIS or 1X PBS) or water can be used for the dilution. Store any thawed aliquot in refrigeration at 2° C to 8° C for up to four weeks, and any frozen aliquot at -20° C to -80° C for up to one year. It is recommended that frozen aliquots be given an amount of standard cryopreservative (such as Ethylene Glycol or Glycerol 5-20% v/v), and refrigerated samples be given an amount of carrier protein (such as heat inactivated FBS or BSA to 0.1% v/v) or non-ionic detergent (such as Triton X-100 or Tween 20 to 0.005% v/v), to aid stability during storage.

SELECT PRODUCT CITATIONS

1. Park, J., et al. 2011. Modulation of CD4⁺ T lymphocyte lineage outcomes with targeted, nanoparticle-mediated cytokine delivery. *Mol. Pharm.* 8: 143-152.
2. Rathje, M., et al. 2011. A peptide derived from the CD loop-D helix region of ciliary neurotrophic factor (CNTF) induces neuronal differentiation and survival by binding to the leukemia inhibitory factor (LIF) receptor and common cytokine receptor chain gp130. *Eur. J. Cell Biol.* 90: 990-999.
3. Lebonvallet, N., et al. 2012. Characterization of neurons from adult human skin-derived precursors in serum-free medium: a PCR array and immunocytochemical analysis. *Exp. Dermatol.* 21: 195-200.
4. Ding, X., et al. 2016. Human GV oocytes generated by mitotically active germ cells obtained from follicular aspirates. *Sci. Rep.* 6: 28218.
5. Salerno-Goncalves, R., et al. 2016. Development of a multicellular three-dimensional organotypic model of the human intestinal mucosa grown under microgravity. *J. Vis. Exp.* 113: 54148.
6. Bataille, A., et al. 2020. *In vitro* differentiation of human skin-derived cells into functional sensory neurons-like. *Cells* 9: 1000.
7. Eghbalsaid, S., et al. 2020. A versatile bulk electrotransfection protocol for murine embryonic fibroblasts and iPS cells. *Sci. Rep.* 10: 13332.
8. Yu, X., et al. 2023. Oocyte arrested at metaphase II stage were derived from human pluripotent stem cells *in vitro*. *Stem Cell Rev. Rep.* 19: 1067-1081.
9. Eghbalsaid, S. and Kues, W.A. 2023. CRISPR/Cas9-mediated targeted knock-in of large constructs using nocodazole and RNase HII. *Sci. Rep.* 13: 2690.

RESEARCH USE

For research use only, not for use in diagnostic procedures.