

DOCK 180 siRNA (h2): sc-44294

BACKGROUND

The v-Crk oncogene product shares homologous amino acid sequences, designated Src homology region 2 (SH2) and SH3, with many molecules involved in signal transduction. The v-Crk cellular homolog, c-Crk, is a member of a newly emerging class of genes including Nck and GRB2/ASH which encode proteins that consist primarily of SH2 and SH3 domains. Two distinct human c-Crk cDNAs, designated Crk I and Crk II, have been identified and shown to represent alternative splice products of c-Crk. The major translational product of c-Crk I has been identified as a variably expressed protein, while c-Crk II encodes a widely expressed protein and a more variably expressed protein. The major c-Crk transforming activity appears associated with c-Crk I p28 expression. DOCK 180, a protein downstream of Crk, has been identified as a major Crk-associated protein. When DOCK 180 is recruited to the plasma membrane from a cytoplasmic reservoir, presumably by Crk, changes in cellular morphology and spindle formation occur, suggesting DOCK 180 to be a Crk effector molecule.

REFERENCES

1. Mayer, B.J., et al. 1990. Association of the v-Crk oncogene product with phosphotyrosine-containing proteins and protein kinase activity. *Proc. Natl. Acad. Sci. USA* 87: 2638-2642.
2. Matsuda, M., et al. 1990. Binding of transforming protein, P47gag-Crk, to a broad range of phosphotyrosine-containing proteins. *Science* 248: 1537-1539.
3. Mayer, B.J., et al. 1990. Mutagenic analysis of the v-Crk oncogene: requirement for SH2 and SH3 domains, and correlation between increased cellular phosphotyrosine and transformation. *J. Virol.* 64: 3581-3589.

CHROMOSOMAL LOCATION

Genetic locus: DOCK1 (human) mapping to 10q26.2.

PRODUCT

DOCK 180 siRNA (h2) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see DOCK 180 shRNA Plasmid (h2): sc-44294-SH and DOCK 180 shRNA (h2) Lentiviral Particles: sc-44294-V as alternate gene silencing products.

For independent verification of DOCK 180 (h2) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-44294A, sc-44294B and sc-44294C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

DOCK 180 siRNA (h2) is recommended for the inhibition of DOCK 180 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

DOCK 180 (E-2): sc-514080 is recommended as a control antibody for monitoring of DOCK 180 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor DOCK 180 gene expression knockdown using RT-PCR Primer: DOCK 180 (h2)-PR: sc-44294-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

1. Li, W., et al. 2018. HGF-induced formation of the Met-Axl-Elmo2-DOCK 180 complex promotes Rac 1 activation, receptor clustering, and cancer cell migration and invasion. *J. Biol. Chem.* 293: 15397-15418.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.