

mucolipin 1 siRNA (m): sc-44520

BACKGROUND

The gene encoding human mucolipin 1 maps to chromosome 19p13.2. Mutations in this gene cause a rare autosomal recessive lysosomal storage disease known as mucopolipidosis type IV (MLIV). Clinical characteristics of MLIV include psychomotor retardation, retinal degeneration, corneal opacities and strabismus. Mucolipin 1 localizes to the plasma membrane and contains six transmembrane domains. The carboxy-terminus of mucolipin 1 shares sequence homology with Polycystin-2 and the transient receptor potential cation channel family. The concentration of intracellular Ca^{2+} regulates the permeability of mucolipin 1 to Ca^{2+} , Na^{+} and K^{+} . The influence of Ca^{2+} on mucolipin 1 represents a possible role for mucolipin 1 in lysosomal exocytosis and the trafficking of late endosomes and lysosomes.

REFERENCES

- Merin, S., Livni, N., Berman, E.R. and Yatziv, S. 1975. Mucopolipidosis IV: ocular, systemic and ultrastructural findings. *Invest. Ophthalmol.* 14: 437-448.
- Bargal, R., Avidan, N., Ben-Asher, E., Olender, Z., Zeigler, M., Frumkin, A., Raas-Rothschild, A., Glusman, G., Lancet, D. and Bach, G. 2000. Identification of the gene causing mucopolipidosis type IV. *Nat. Genet.* 26: 118-123.
- Sun, M., Goldin, E., Stahl, S., Falardeau, J.L., Kennedy, J.C., Bove, C., Acierno, J.S., Jr., Kaneski, C.R., Nagle, J., Bromley, M.C., Colman, M., Schiffmann, R. and Slaugenhaupt, S.A. 2000. Mucopolipidosis type IV is caused by mutations in a gene encoding a novel transient receptor potential channel. *Hum. Mol. Genet.* 9: 2471-2478.
- Bassi, M.T., Manzoni, M., Monti, E., Pizzo, M.T., Ballabio, A. and Borsani, G. 2000. Cloning of the gene encoding a novel integral membrane protein, mucolipidin—and identification of the two major founder mutations causing mucopolipidosis type IV. *Am. J. Hum. Genet.* 67: 1110-1120.
- LaPlante, J.M., Falardeau, J., Sun, M., Kanazirska, M., Brown, E.M., Slaugenhaupt, S.A. and Vassilev, P.M. 2002. Identification and characterization of the single channel function of human mucolipin-1 implicated in mucopolipidosis type IV, a disorder affecting the lysosomal pathway. *FEBS Lett.* 532: 183-187.

CHROMOSOMAL LOCATION

Genetic locus: Mcoln1 (mouse) mapping to 8 A1.1.

PRODUCT

mucolipin 1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see mucolipin 1 shRNA Plasmid (m): sc-44520-SH and mucolipin 1 shRNA (m) Lentiviral Particles: sc-44520-V as alternate gene silencing products.

For independent verification of mucolipin 1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-44520A, sc-44520B and sc-44520C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

mucolipin 1 siRNA (m) is recommended for the inhibition of mucolipin 1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor mucolipin 1 gene expression knockdown using RT-PCR Primer: mucolipin 1 (m)-PR: sc-44520-PR (20 μ l, 538 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

- Park, H.J., Gholam-Zadeh, M., Suh, J.H. and Choi, H.S. 2019. Lycorine attenuates autophagy in osteoclasts via an axis of mROS/TRPML1/TFEB to reduce LPS-induced bone loss. *Oxid. Med. Cell. Longev.* 2019: 8982147.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.