

MYL3 (MLM508): sc-47719

BACKGROUND

Myosin, the major component of thick muscle filaments, is a long asymmetric molecule containing a globular head and a long tail. Activation of smooth and cardiac/ventricular muscle primarily involves pathways which increase calcium and myosin phosphorylation, resulting in contraction. Myosin in vertebrate striated muscle is composed of two heavy chains and four light chains. There are two distinct types of light chains: the phosphorylatable, regulatory or MLC2 type; and the nonphosphorylatable, alkali or MLC1 and MLC3 types. Myosin light chain phosphatase acts to regulate muscle contraction by dephosphorylating activated Myosin light chain. The role of Myosin alkali light chains in vertebrate skeletal muscle is poorly understood, although alkali light chains in smooth muscle may be involved with the active site of Myosin. Several isoforms of Myosin alkali light chains have been identified, encoded by a family of Myosin light chain genes. Each is associated with different muscle types. Human Myosin light chain can be used as a cardiac marker. Myosin light chain 3, encoded by MYL3, is an alkali light chain also referred to as both the ventricular isoform (MLC1v) and slow skeletal muscle isoform. Myosin light chain 3 proteins in human and mouse share 91% sequence identity overall.

REFERENCES

1. Barton, P.J. and Buckingham, M.E. 1985. The Myosin alkali light chain proteins and their genes. *Biochem. J.* 231: 249-261.
2. Seidel, U., et al. 1987. The complete nucleotide sequences of cDNA clones coding for human Myosin light chains 1 and 3. *Nucleic Acids Res.* 15: 4989.
3. Cohen-Hagueneauer, O., et al. 1988. Assignment of the human fast skeletal muscle Myosin alkali light chains gene (MLC1F/MLC3F) to 2q32.1-2qter. *Hum. Genet.* 78: 65-70.
4. Katoh, H., et al. 1992. Development of an immunoradiometric assay kit for ventricular Myosin light chain I with monoclonal antibodies. *Clin. Chem.* 38: 170-171.
5. Sanbe, A., et al. 1999. Abnormal cardiac structure and function in mice expressing nonphosphorylatable cardiac regulatory Myosin light chain 2. *J. Biol. Chem.* 274: 21085-21094.
6. Davis, J.S., et al. 2001. The overall pattern of cardiac contraction depends on a spatial gradient of Myosin regulatory light chain phosphorylation. *Cell* 107: 631-641.
7. Yamashita, H., et al. 2003. Myosin light chain isoforms modify force-generating ability of cardiac Myosin by changing the kinetics of Actin-Myosin interaction. *Cardiovasc. Res.* 60: 580-588.
8. Bicer, S. and Reiser, P.J. 2004. Myosin light chain 1 isoforms in slow fibers from global and orbital layers of canine rectus muscles. *Invest. Ophthalmol. Vis. Sci.* 45: 138-143.

CHROMOSOMAL LOCATION

Genetic locus: MYL3 (human) mapping to 3p21.31.

SOURCE

MYL3 (MLM508) is a mouse monoclonal antibody raised against purified cardiac MYL3.

PRODUCT

Each vial contains 100 µg IgG_{2a} in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

MYL3 (MLM508) is recommended for detection of myosin light chain 3 of human and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500).


MYL3 (MLM508) is also recommended for detection of myosin light chain 3 in additional species, including bovine and porcine.

Suitable for use as control antibody for MYL3 siRNA (h): sc-44542, MYL3 shRNA Plasmid (h): sc-44542-SH and MYL3 shRNA (h) Lentiviral Particles: sc-44542-V.


Molecular Weight of MYL3: 25 kDa.

Positive Controls: rat heart extract: sc-2393.

DATA



MYL3 (MLM508): sc-47719. Western blot analysis of myosin light chain expression in rat heart tissue extract.



MYL3 (MLM508): sc-47719. Immunoperoxidase staining of formalin fixed, paraffin-embedded human heart muscle tissue showing cytoplasmic staining of myocytes at low (A) and high (B) magnification. Kindly provided by The Swedish Human Protein Atlas (HPA) program.

SELECT PRODUCT CITATIONS

1. Lindskog, C., et al. 2015. The human cardiac and skeletal muscle proteomes defined by transcriptomics and antibody-based profiling. *BMC Genomics* 16: 475.
2. Gao, R., et al. 2019. Pioneering function of Isl1 in the epigenetic control of cardiomyocyte cell fate. *Cell Res.* 29: 486-501.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.