

RANTES (mBA-68): sc-4952

BACKGROUND

Structurally, C-C or β -chemokines are characterized by a set of conserved, adjacent cysteines. Members of this family include MCP-1, MCP-2, MCP-3, MIP-1 α , MIP-1 β , RANTES and I-309. RANTES (regulated upon activation, normal T cell expressed and secreted) is expressed by platelets, eosinophils, fibroblasts, macrophages, endothelial cells and T lymphocytes. Consistent with its belonging to the chemokine family, RANTES exhibits strong chemoattractant activity towards monocytes and NK cells. I-309 was initially identified as a factor present in γ/δ T lymphocytes. I-309 cDNA encodes a protein 73 amino acids in length with one potential N-linked glycosylation site and migrates at a molecular weight of 15 kDa. Unlike the other members of the C-C family, I-309 does not induce chemotaxis in natural killer (NK) cells.

REFERENCES

1. Miller, M.D., Hata, S., De Waal Malefyt, R., and Krangel, M.S. 1989. A novel polypeptide secreted by activated human T lymphocytes. *J. Immunol.* 143: 2907-2916.
2. Lloyd, A.R., Oppenheim, J.J., Kelvin, D.J., and Taub, D.D. 1996. Chemokines regulate T cell adherence to recombinant adhesion molecules and extracellular matrix proteins. *J. Immunol.* 156: 932-938.
3. Taub, D.D., Ortaldo, J.R., Turcovski-Corrales, S.M., Key, M.L., Longo, D.L., and Murphy, W.J. 1996. Beta chemokines costimulate lymphocyte cytolysis, proliferation and lymphokine production. *J. Leukoc. Biol.* 59: 81-89.
4. Wells, T.N., Power, C.A., Lusti-Narasimhan, M., Hoogewerf, A.J., Cooke, R.M., Chung, C.W., Peitsch, M.C., and Proudfoot, A.E. 1996. Selectivity and antagonism of chemokine receptors. *J. Leukoc. Biol.* 59: 53-60.
5. Wang, J.H., Devalia, J.L., Xia, C., Sapsford, R.J., and Davies, R.J. 1996. Expression of RANTES by human bronchial epithelial cells *in vitro* and *in vivo* and the effect of corticosteroids. *Am. J. Respir. Cell Mol. Biol.* 14: 27-35.
6. Ying, S., Meng, Q., Taborda-Barata, L., Corrigan, C.J., Barkans, J., Assoufi, B., Moqbel, R., Durham, S.R., and Kay, A.B. 1996. Human eosinophils express messenger RNA encoding RANTES and store and release biologically active RANTES protein. *Euro. J. Immunol.* 26: 70-76.
7. Loetscher, P., Seitz, M., Clark-Lewis, I., Baggiolini, M., and Moser, B. 1996. Activation of NK cells by C-C chemokines. Chemotaxis, Ca²⁺ mobilization and enzyme release. *J. Immunol.* 156: 322-327.

SOURCE

RANTES (mBA-68) is produced in *E. coli* as 35 kDa biologically active, tagged fusion protein corresponding to 68 amino acids of RANTES of mouse origin.

PRODUCT

RANTES (mBA-68) is purified from bacterial lysates (>98%); supplied as 50 μ g purified protein.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

BIOLOGICAL ACTIVITY

RANTES (mBA-68) is biologically active as determined by its ability to chemoattract total human lymphocyte population and total murine T cell population using a concentration range of 1.0–10 mg/ml.

RECONSTITUTION

In order to avoid freeze/thaw damaging of the active protein, dilute protein when first used to desired working concentration. Either a sterile filtered standard buffer (such as 50mM TRIS or 1X PBS) or water can be used for the dilution. Store any thawed aliquot in refrigeration at 2° C to 8° C for up to four weeks, and any frozen aliquot at -20° C to -80° C for up to one year. It is recommended that frozen aliquots be given an amount of standard cryopreservative (such as Ethylene Glycol or Glycerol 5-20% v/v), and refrigerated samples be given an amount of carrier protein (such as heat inactivated FBS or BSA to 0.1% v/v) or non-ionic detergent (such as Triton X-100 or Tween 20 to 0.005% v/v), to aid stability during storage.

STORAGE

Store desiccated at -20° C; stable for one year from the date of shipment.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.