

LARS2 (D-10): sc-514756

BACKGROUND

LARS2 (leucyl-tRNA synthetase 2, mitochondrial) is also known as LEURS (leucine-tRNA ligase) and is a 903 amino acid protein. LARS2 is a member of the class-I aminoacyl-tRNA synthetase family and is localized to the mitochondrial matrix. LARS2 catalyzes the aminoacylation of leucine to tRNA^{Leu} via a two step reaction during protein synthesis. The two step reaction begins by LARS2 activating leucine with an ATP molecule which yields an adenylate intermediate that then transfers the activated leucine to the 3'-end of the target tRNA. tRNA^{Leu} has a variable loop with a specific sequence and orientation which is thought to be important for interaction with LARS2. LARS2 is upregulated in bipolar disorder and schizophrenia and is thought to be over-expressed in an attempt to cause a mutated tRNA^{Leu}, tRNA^{Leu(UUR)}, to go through aminoacylation. Diabetes is also thought to be associated with upregulation of LARS2 which may promote intolerance of glucose.

REFERENCES

1. Han, W., et al. 2001. Gene cloning, expression and purification of human mitochondrial tRNA^{Leu(UUR)} and its mutant. *Sci. China, C, Life Sci.* 44: 113-120.
2. Munakata, K., et al. 2005. Mitochondrial DNA 3243A>G mutation and increased expression of LARS2 gene in the brains of patients with bipolar disorder and schizophrenia. *Biol. Psychiatry* 57: 525-532.
3. 't Hart, L.M., et al. 2005. Evidence that the mitochondrial leucyl tRNA synthetase (LARS2) gene represents a novel type 2 diabetes susceptibility gene. *Diabetes* 54: 1892-1895.
4. Zhai, Y., et al. 2007. Modulation of substrate specificity within the amino acid editing site of leucyl-tRNA synthetase. *Biochemistry* 46: 3331-3337.
5. Lue, S.W. and Kelley, S.O. 2007. A single residue in leucyl-tRNA synthetase affecting amino acid specificity and tRNA aminoacylation. *Biochemistry* 46: 4466-4472.
6. Fukunaga, R. and Yokoyama, S. 2007. The C-terminal domain of the archaeal leucyl-tRNA synthetase prevents misediting of isoleucyl-tRNA^{Ile}. *Biochemistry* 46: 4985-4996.

CHROMOSOMAL LOCATION

Genetic locus: LARS2 (human) mapping to 3p21.31.

SOURCE

LARS2 (D-10) is a mouse monoclonal antibody specific for an epitope mapping between amino acids 217-232 within an internal region of LARS2 of human origin.

PRODUCT

Each vial contains 200 µg IgM kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-514756 P, (100 µg peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% stabilizer protein).

APPLICATIONS

LARS2 (D-10) is recommended for detection of LARS2 of human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for LARS2 siRNA (h): sc-78462, LARS2 shRNA Plasmid (h): sc-78462-SH and LARS2 shRNA (h) Lentiviral Particles: sc-78462-V.


Molecular Weight of LARS2: 102 kDa.

Positive Controls: Raji whole cell lysate: sc-364236, CCRF-CEM cell lysate: sc-2225 or Hep G2 cell lysate: sc-2227.


RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein L-Agarose: sc-2336 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

DATA



LARS2 (D-10): sc-514756. Western blot analysis of LARS2 expression in K-562 (A), Jurkat (B), Raji (C), MOLT-4 (D), Hep G2 (E) and CCRF-CEM (F) whole cell lysates.



LARS2 (D-10): sc-514756. Western blot analysis of LARS2 expression in HeLa whole cell lysate.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.