

DYNLL1 (AT13C9): sc-517399

BACKGROUND

Dyneins are multisubunit, high molecular weight ATPases that interact with microtubules to generate force by converting the chemical energy of ATP into the mechanical energy of movement. Cytoplasmic or axonemal Dynein heavy, intermediate, light and light-intermediate chains are all components of minus end-directed motors; the complex transports cellular cargos towards the central region of the cell. The highly conserved DYNLL proteins were originally identified as light chains for microtubule-based motor protein Dynein. In mammals there are two closely related isoforms expressed, DYNLL1 and DYNLL2 which share 93% sequence identity at the protein level. DYNLL1 (Dynein light chain 1) also designated, DLC8 or PIN (Protein inhibitor of neuronal nitric oxide synthase) has been identified as a protein that interacts with NOS1 resulting in NOS1 inhibition. Dimerization is required for NOS1 activity and DYNLL1 has been shown to destabilize the NOS1 dimer. Nitric oxide may be involved in several processes such as apoptosis, synaptogenesis and neuronal development; thus DYNLL1 is implicated in these processes as well. DYNLL1 is a ubiquitously expressed protein that exhibits high expression in testis and moderate expression in brain. DYNLL2 (Dynein light chain 2) is subject to a unique alternative splicing event which is implicated in Myosin Va binding specificity.

REFERENCES

1. Dick, T., Ray, K., Salz, H.K. and Chia, W. 1996. Cytoplasmic dynein (ddlc1) mutations cause morphogenetic defects and apoptotic cell death in *Drosophila melanogaster*. Mol. Cell. Biol. 16: 1966-1977.
2. Jaffrey, S.R. and Snyder, S.H. 1996. PIN: an associated protein inhibitor of neuronal nitric oxide synthase. Science 274: 774-777.
3. Fuhrmann, J.C., Kins, S., Rostaing, P., El Far, O., Kirsch, J., Sheng, M., Triller, A., Betz, H. and Kneussel, M. 2002. Gephyrin interacts with Dynein light chains 1 and 2, components of motor protein complexes. J. Neurosci. 22: 5393-5402.
4. Asai, D.J. and Wilkes, D.E. 2004. The dynein heavy chain family. J. Eukaryot. Microbiol. 51: 23-29.
5. Mallik, R., Carter, B.C., Lex, S.A., King, S.J. and Gross, S.P. 2004. Cytoplasmic dynein functions as a gear in response to load. Nature 427: 649-652.
6. Lee, W.L., Kaiser, M.A. and Cooper, J.A. 2005. The offloading model for dynein function: differential function of motor subunits. J. Cell Biol. 168: 201-207.
7. Wagner, W., Fodor, E., Ginsburg, A. and Hammer, J.A. 2006. The binding of DYNLL2 to myosin Va requires alternatively spliced exon B and stabilizes a portion of the myosin's coiled-coil domain. Biochemistry 45: 11564-11577.
8. Pfister, K.K., Shah, P.R., Hummerich, H., Russ, A., Cotton, J., Annuar, A.A., King, S.M. and Fisher, E.M. 2006. Genetic analysis of the cytoplasmic dynein subunit families. PLoS Genet. 2: e1.

CHROMOSOMAL LOCATION

Genetic locus: DYNLL1 (human) mapping to 12q24.31.

SOURCE

DYNLL1 (AT13C9) is a mouse monoclonal antibody raised against a recombinant protein corresponding to amino acids 1-89 of DYNLL1 of human origin.

PRODUCT

Each vial contains 100 µg IgG_{2b} kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

DYNLL1 (AT13C9) is recommended for detection of DYNLL1 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), flow cytometry (1 µg per 1 x 10⁶ cells) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for DYNLL1 siRNA (h): sc-36228, DYNLL1 shRNA Plasmid (h): sc-36228-SH and DYNLL1 shRNA (h) Lentiviral Particles: sc-36228-V.

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.