

Rab 13 (160A2F): sc-517609

BACKGROUND

The Ras-related superfamily of guanine nucleotide binding proteins, which includes the R-Ras, Rap, Ral/Rec and Rho/Rab superfamilies, exhibits 30-60% homology with Ras p21. Accumulating data suggests an important role for Rab proteins, either in endocytosis or in biosynthetic protein transport. The Rab family of small G proteins play an important role in determining the specificity of vesicular transport pathways. Rab 13 and Rab 3B localize to tight junctions in epithelial cells and to cytoplasmic vesicular structures in cells lacking tight junctions. Rab 13 can be detected in the junctional complex regions of a variety of epithelia, including intestine, kidney and liver.

REFERENCES

1. Novick, P. and Brennwald, P. 1993. Friends and family: the role of the Rab GTPases in vesicular traffic. *Cell* 75: 597-601.
2. Ferro-Novick, S. and Novick, P. 1993. The role of GTP-binding proteins in transport along the exocytic pathway. *Annu. Rev. Cell Biol.* 9: 575-599.
3. Zahraoui, A., Joberty, G., Arpin, M., Fontaine, J.J., Hellio, R., Tavitian, A. and Louvard, D. 1994. A small Rab GTPase is distributed in cytoplasmic vesicles in non polarized cells but colocalizes with the tight junction marker ZO-1 in polarized epithelial cells. *J. Cell Biol.* 124: 101-115.
4. Leek, J.P., Hamlin, P.J., Wilton, J. and Lench, N.J. 1997. Assignment of the Rab13 gene (RAB13) to human chromosome band 12q13 by *in situ* hybridization. *Cytogenet. Cell Genet.* 79: 210-211.
5. Sheth, B., Fontaine, J.J., Ponza, E., McCallum, A., Page, A., Citi, S., Louvard, D., Zahraoui, A. and Fleming, T.P. 2000. Differentiation of the epithelial apical junctional complex during mouse preimplantation development: a role for Rab13 in the early maturation of the tight junction. *Mech. Dev.* 97: 93-104.
6. Yamamoto, Y., Nishimura, N., Morimoto, S., Kitamura, H., Manabe, S., Kanayama, H.O., Kagawa, S. and Sasaki, T. 2003. Distinct roles of Rab 3B and Rab 13 in the polarized transport of apical, basolateral and tight junctional membrane proteins to the plasma membrane. *Biochem. Biophys. Res. Commun.* 308: 270-275.
7. Morimoto, S., Nishimura, N., Terai, T., Manabe, S., Yamamoto, Y., Shinahara, W., Miyake, H., Tashiro, S., Shimada, M. and Sasaki T. 2005. Rab13 mediates the continuous endocytic recycling of occludin to the cell surface. *J. Biol. Chem.* 280:2220-2228.
8. Marzesco, A.M. and Zahraoui, A. 2005. Assay of Rab13 in regulating epithelial tight junction assembly. *Meth. Enzymol.* 403:182-193.

CHROMOSOMAL LOCATION

Genetic locus: RAB13 (human) mapping to 1q21.3.

SOURCE

Rab 13 (160A2F) is a mouse monoclonal antibody raised against amino acids 1-203 representing full length Rab 13 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 100 µg IgG₁ in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

Rab 13 (160A2F) is recommended for detection of Rab 13 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000).

Suitable for use as control antibody for Rab 13 siRNA (h): sc-44060, Rab 13 shRNA Plasmid (h): sc-44060-SH and Rab 13 shRNA (h) Lentiviral Particles: sc-44060-V.

Molecular Weight of Rab 13: 25 kDa.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.