# 14-3-3 ζ (G-2): sc-518031



The Power to Question

#### **BACKGROUND**

14-3-3 proteins regulate many cellular processes relevant to cancer biology, notably apoptosis, mitogenic signaling and cell-cycle checkpoints. Seven isoforms comprise this family of signaling intermediates, denoted 14-3-3  $\beta,\gamma,\epsilon,\zeta,\eta,\theta$  and  $\sigma$ . 14-3-3 proteins form dimers that present two binding sites for ligand proteins, thereby bringing together two proteins that may not otherwise associate. These ligands largely share a 14-3-3 consensus binding motif and exhibit serine/threonine phosphorylation. 14-3-3 proteins function in broad regulation of these ligand proteins, by cytoplasmic sequestration, occupation of interaction domains and import/export sequences, prevention of degradation, activation/repression of enzymatic activity and facilitation of protein modification, and thus loss of expression contributes to a vast array of pathogenic cellular activities.

#### **CHROMOSOMAL LOCATION**

Genetic locus: YWHAZ (human) mapping to 8q22.3, YWHAE (human) mapping to 17p13.3; Ywhaz (mouse) mapping to 15 B3.1, Ywhae (mouse) mapping to 11 B5.

#### **SOURCE**

14-3-3  $\zeta$  (G-2) is a mouse monoclonal antibody specific for an epitope mapping between amino acids 122-140 within an internal region of 14-3-3  $\zeta$  of human origin.

## **PRODUCT**

Each vial contains 200  $\mu g$   $lgG_3$  kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

14-3-3 ζ (G-2) is available conjugated to agarose (sc-518031 AC), 500 μg/ 0.25 ml agarose in 1 ml, for IP; to HRP (sc-518031 HRP), 200 μg/ml, for WB, IHC(P) and ELISA; to either phycoerythrin (sc-518031 PE), fluorescein (sc-518031 FITC), Alexa Fluor® 488 (sc-518031 AF488), Alexa Fluor® 546 (sc-518031 AF546), Alexa Fluor® 594 (sc-518031 AF594) or Alexa Fluor® 647 (sc-518031 AF647), 200 μg/ml, for WB (RGB), IF, IHC(P) and FCM; and to either Alexa Fluor® 680 (sc-518031 AF680) or Alexa Fluor® 790 (sc-518031 AF790), 200 μg/ml, for Near-Infrared (NIR) WB, IF and FCM.

#### **STORAGE**

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### **APPLICATIONS**

14-3-3  $\zeta$  (G-2) is recommended for detection of 14-3-3  $\zeta$  and, to a lesser extent, 14-3-3  $\epsilon$  of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2  $\mu$ g per 100-500  $\mu$ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


Molecular Weight of 14-3-3 ζ: 30 kDa.


Positive Controls: TK-1 whole cell lysate: sc-364798, M1 whole cell lysate: sc-364782 or Neuro-2A whole cell lysate: sc-364185.

#### **RECOMMENDED SUPPORT REAGENTS**

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-lgGκ BP-HRP: sc-516102 or m-lgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz\* Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-lgGκ BP-FITC: sc-516140 or m-lgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz\* Mounting Medium: sc-24941 or UltraCruz\* Hard-set Mounting Medium: sc-359850.

#### **DATA**





14-3-3 \( (G-2)\) HRP: sc-518031\) HRP. Direct western blot analysis of 14-3-3 \( \) expression in TK-1 (A), M1 (B), Sol8 (C), HL-60 (D), NIH/3T3 (E) and SK-MEL-28 (F) whole cell Ivsates.

14-3-3  $\zeta$  (G-2): sc-518031. Western blot analysis of 14-3-3  $\zeta$  expression in Neuro-2A (**A**), BC $_3$ H1 (**B**), EOC 20 (**C**) and JC (**D**) whole cell lysates.

#### **SELECT PRODUCT CITATIONS**

- Abdrabou, A., et al. 2020. Differential subcellular distribution and translocation of seven 14-3-3 isoforms in response to EGF and during the cell cycle. Int. J. Mol. Sci. 21: 318.
- 2. Mei, J., et al. 2021. YWHAZ interacts with DAAM1 to promote cell migration in breast cancer. Cell Death Discov. 7: 221.
- 3. Huo, Q., et al. 2022. Clinicopathological features and prognostic evaluation of UBR5 in liver cancer patients. Pathol. Oncol. Res. 28: 1610396.

### **RESEARCH USE**

For research use only, not for use in diagnostic procedures.

#### **PROTOCOLS**

See our web site at www.scbt.com for detailed protocols and support products.

Alexa Fluor® is a trademark of Molecular Probes, Inc., Oregon, USA