

Tim10 (C-7): sc-518241

BACKGROUND

The majority of mitochondrial-directed proteins are encoded by the nuclear genome and are transported to the mitochondria via regulated processes involving the mitochondrial Tom and Tim proteins. The mitochondrial Tim protein family is comprised of a large group of evolutionarily conserved proteins that are found in most eukaryotes. Import of nuclear-encoded precursor proteins into and across the mitochondrial inner membrane is mediated by two distinct complexes, the Tim23 complex and the Tim22 complex, which differ in their substrate specificity. Defects in Tim proteins are implicated in several neuro-degenerative diseases, suggesting important roles for Tim proteins in development and health. Tim10, which maps to human chromosome 11q12.1-q12.3, forms heteromeric complexes with Tim9 and Tim12. One complex contains Tim9 and Tim10, which cross-links to the carboxy-terminal domain of Tim23. The carboxy-terminal domain of Tim23 carries all the targeting signals for Tim23, suggesting important role for the Tim9-Tim10 complex in Tim23 import. The other complex contains Tim9, Tim10 and Tim12, which associates with Tim22.

REFERENCES

1. Jin, H., et al. 1999. The human family of Deafness/Dystonia peptide (DDP) related mitochondrial import proteins. *Genomics* 61: 259-267.
2. Bauer, M.F., et al. 1999. The mitochondrial TIM22 preprotein translocase is highly conserved throughout the eukaryotic kingdom. *FEBS Lett.* 464: 41-47.
3. Rassow, J., et al. 1999. The preprotein translocase of the mitochondrial inner membrane: function and evolution. *J. Mol. Biol.* 286: 105-120.
4. Adam, A., et al. 1999. Tim9, a new component of the TIM22.54 translocase in mitochondria. *EMBO J.* 18: 313-319.
5. Davis, A.J., et al. 2000. Two intermembrane space TIM complexes interact with different domains of Tim23p during its import into mitochondria. *J. Cell Biol.* 150: 1271-1282.
6. Bauer, M.F. and Neupert, W. 2001. Import of proteins into mitochondria: a novel pathomechanism for progressive neurodegeneration. *J. Inher. Metab. Dis.* 24: 166-180.
7. LocusLink Report (LocusID: 26519). <http://www.ncbi.nlm.nih.gov/LocusLink/>

CHROMOSOMAL LOCATION

Genetic locus: TIMM10 (human) mapping to 11q12.1.

SOURCE

Tim10 (C-7) is a mouse monoclonal antibody specific for an epitope mapping between amino acids 52-71 of Tim10 of human origin.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Each vial contains 200 µg IgG₁ kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Tim10 (C-7) is available conjugated to agarose (sc-518241 AC), 500 µg/0.25 ml agarose in 1 ml, for IP; to HRP (sc-518241 HRP), 200 µg/ml, for WB, IHC(P) and ELISA; to either phycoerythrin (sc-518241 PE), fluorescein (sc-518241 FITC), Alexa Fluor® 488 (sc-518241 AF488), Alexa Fluor® 546 (sc-518241 AF546), Alexa Fluor® 594 (sc-518241 AF594) or Alexa Fluor® 647 (sc-518241 AF647), 200 µg/ml, for WB (RGB), IF, IHC(P) and FCM; and to either Alexa Fluor® 680 (sc-518241 AF680) or Alexa Fluor® 790 (sc-518241 AF790), 200 µg/ml, for Near-Infrared (NIR) WB, IF and FCM.

Alexa Fluor® is a trademark of Molecular Probes, Inc., Oregon, USA

APPLICATIONS


Tim10 (C-7) is recommended for detection of Tim10 of human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for Tim10 siRNA (h): sc-41255, Tim10 shRNA Plasmid (h): sc-41255-SH and Tim10 shRNA (h) Lentiviral Particles: sc-41255-V.

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

DATA



Tim10 (C-7): sc-518241. Western blot analysis of human recombinant Tim10 fusion protein. Detection reagent used: m-IgG₁ BP-HRP: sc-525408.

RESEARCH USE

For research use only, not for use in diagnostic procedures.