

COL4A3 (1042): sc-52317

BACKGROUND

The extensive family of COL gene products (collagens) is composed of several chain types, including fibril-forming interstitial Collagens (Types I, II, III and V) and basement membrane Collagens (Type IV), each type containing multiple isoforms. Collagens are fibrous, extracellular matrix proteins with high tensile strength and are the major components of connective tissue, such as tendons and cartilage. All collagens contain a triple helix domain and frequently show lateral self-association in order to form complex connective tissues. Several collagens also play a role in cell adhesion, important for maintaining normal tissue architecture and function.

REFERENCES

- Bateman, J.F., et al. 1996. Collagen Superfamily. In Comper, W.D., ed., Extracellular Matrix. Amsterdam: Harwood Academic Publishers, 22-67.
- McCarthy, J.B., et al. 1996. Cell adhesion to collagenous matrices. *Biopolymers* 40: 371-381.
- Engel, J. 1997. Versatile collagens in invertebrates. *Science* 277: 1785-1786.
- Cremer, M.A., et al. 1998. The cartilage collagens: a review of their structure, organization and role in the pathogenesis of experimental arthritis in animals and in human rheumatic disease. *J. Mol. Med.* 76: 275-288.
- Boskey, A.L., et al. 1999. Collagen and bone strength. *J. Bone Miner. Res.* 14: 330-335.
- Alberio, L. and Dale, G.L. 1999. Platelet-collagen interactions: membrane receptors and intracellular signaling pathways. *Eur. J. Clin. Invest.* 29: 1066-1076.
- Kalluri, R. 2002. Discovery of type IV Collagen non-collagenous domains as novel integrin ligands and endogenous inhibitors of angiogenesis. *Cold Spring Harb. Symp. Quant. Biol.* 67: 255-266.
- Pescucci, C., et al. 2003. Type-IV Collagen related diseases. *J. Nephrol.* 16: 314-316.
- Kim, S.J., et al. 2004. Effects of type IV Collagen and laminin on the cryopreservation of human embryonic stem cells. *Stem Cells* 22: 950-961.

CHROMOSOMAL LOCATION

Genetic locus: COL4A3 (human) mapping to 2q36.3.

SOURCE

COL4A3 (1042) is a mouse monoclonal antibody raised against Collagen α 3 Type IV from isolated placenta of human origin.

PRODUCT

Each vial contains 500 μ l ascites containing IgG_{2b} with < 0.1% sodium azide.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

APPLICATIONS

COL4A3 (1042) is recommended for detection of Collagen α 3 Type IV of human origin by Western Blotting (starting dilution to be determined by researcher, dilution range 1:10-1:200), immunofluorescence (starting dilution to be determined by researcher, dilution range 1:10-1:200) and immunohistochemistry (including paraffin-embedded sections) (starting dilution to be determined by researcher, dilution range 1:10-1:200).

Molecular Weight of COL4A3: 160-190 kDa.

Positive Controls: Jurkat whole cell lysate: sc-2204, HeLa whole cell lysate: sc-2200 or Hs68 cell lysate: sc-2230.

SELECT PRODUCT CITATIONS

- Fu, J., et al. 2009. Paeoniflorin prevents diabetic nephropathy in rats. *Comp. Med.* 59: 557-566.
- Zhang, W., et al. 2010. Three dimensional dental epithelial-mesenchymal constructs of predetermined size and shape for tooth regeneration. *Biomaterials* 31: 7995-8003.
- Lambers, C., et al. 2013. The interaction of endothelin-1 and TGF- β 1 mediates vascular cell remodeling. *PLoS ONE* 8: e73399.
- Barreto-Ortiz, S.F., et al. 2013. A novel *in vitro* model for microvasculature reveals regulation of circumferential ECM organization by curvature. *PLoS ONE* 8: e81061.
- Liang, D., et al. 2015. EGFR inhibition protects cardiac damage and remodeling through attenuating oxidative stress in STZ-induced diabetic mouse model. *J. Mol. Cell. Cardiol.* 82: 63-74.
- Laperle, A., et al. 2015. α -5 laminin synthesized by human pluripotent stem cells promotes self-renewal. *Stem Cell Reports* 5: 195-206.
- Yang, X., et al. 2021. Differential expression of COL4A3 and collagen in upward and downward progressing types of nasopharyngeal carcinoma. *Oncol. Lett.* 21: 223.
- Barbon, S., et al. 2022. Preclinical development of bioengineered allografts derived from decellularized human diaphragm. *Biomedicines* 10: 739.
- Barbon, S., et al. 2022. Bio-engineered scaffolds derived from decellularized human esophagus for functional organ reconstruction. *Cells* 11: 2945.

STORAGE

For immediate and continuous use, store at 4° C for up to one month. For sporadic use, freeze in working aliquots in order to avoid repeated freeze/thaw cycles. If turbidity is evident upon prolonged storage, clarify solution by centrifugation.

RESEARCH USE

For research use only, not for use in diagnostic procedures.