



m-IgG₁ BP-B: sc-533660

BACKGROUND

Mouse IgG₁ binding protein (IgG₁ BP) conjugated to biotin (B) is a strongly recommended alternative to conventional goat/rabbit anti-mouse IgG secondary antibodies for immunohistochemical (IHC) signal enhancement. Mouse IgG₁ binding protein is a highly specific reagent that provides strong signal with minimal background and virtually complete elimination of lot to lot variation associated with conventionally generated secondary antibodies. Mouse IgG₁ binding protein (m-IgG₁ BP) is suitable for binding to most, but not all, mouse IgG₁ immunoglobulins, comprising approximately 55% of SCBT's mouse monoclonal antibodies; not suitable for use with mouse monoclonal IgG_{2a}, IgG_{2b}, IgG₃, IgM, IgA and IgE antibodies. It may slightly cross react with mouse IgG_{2b} or goat IgG antibodies. Not cross reactive with human, rat or rabbit IgG antibodies.

SOURCE

m-IgG₁ BP-B is a purified recombinant mouse IgG₁ binding protein conjugated to biotin (B).

PRODUCT

Each vial contains 100 µg mouse IgG₁ binding protein-B in 0.5 ml of PBS containing 1% stabilizer protein and 0.02% sodium azide.

APPLICATIONS

m-IgG₁ BP-B is recommended for detection of mouse IgG₁ by immunohistochemistry (starting dilution: 1:50, dilution range 1:50-1:200). Optimal dilution to be determined by titration.

RECOMMENDED SUPPORT PRODUCTS

- CrystalCruz™ Cover Glasses, 22 x 50 mm, precleaned: sc-24975
- PBS (Phosphate Buffered Saline), powder, 1 packet: sc-24947
- Formaldehyde, 37% formaldehyde solution, 25 ml: sc-203049
- Hydrogen Peroxide, 30% solution, 100 ml: sc-203336
- Organo/Limonene Mount, non-toxic alternative to Permount, 100 ml: sc-45087
- UltraCruz® Mounting Medium, aqueous-based, 10 ml: sc-24941
- ImmunoHistoMount, aqueous-based mounting medium, 30 ml: sc-45086
- Immuno In Situ Mount, for use with *in situ* hybridization, 30 ml: sc-45088
- Paraffin, for the preparation of tissue samples for staining, 500 g: sc-286633
- Xylenes, mixed isomers with ethylbenzene, 500 ml: sc-237422
- Hematoxylin, Gill's Formulation #2; nuclear counter stain, 100 ml: sc-24973

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.