



m-IgG_{2b} BP-FITC: sc-542743

BACKGROUND

Mouse IgG_{2b} binding protein (m-IgG_{2b} BP) conjugated to fluorescein isothiocyanate (FITC) is a strongly recommended alternative to conventional goat/rabbit anti-mouse IgG secondary antibodies for immunofluorescence (IF) and flow cytometry (FCM) signal enhancement. Mouse IgG_{2b} binding protein is a highly specific reagent that provides strong signal with minimal background and virtually complete elimination of lot to lot variation associated with conventionally generated secondary antibodies. Mouse IgG_{2b} binding protein (m-IgG_{2b} BP) is suitable for binding to most, but not all mouse monoclonal IgG_{2b} antibodies; not suitable for use with mouse monoclonal IgG₁, IgG_{2a}, IgG₃, IgM, IgA or IgE antibodies. Not cross reactive with human, rat or goat IgG antibodies.

SOURCE

m-IgG_{2b} BP-FITC is a purified recombinant mouse IgG_{2b} binding protein conjugated to fluorescein isothiocyanate (FITC).

PRODUCT

Each vial contains 50 µg mouse IgG_{2b} binding protein-FITC in 0.5 ml of PBS containing 0.1% gelatin and 0.1% sodium azide.

APPLICATIONS

m-IgG_{2b} BP-FITC is recommended for detection of mouse IgG_{2b} by immunofluorescence (starting dilution 1:50, dilution range 1:50-1:200) and flow cytometry (0.5-1 µg per 1 x 10⁶ cells). Optimal dilution to be determined by titration.

FITC excitation/emission spectrum peak wavelengths: 495 nm/519 nm

RECOMMENDED SUPPORT PRODUCTS

- CrystalCruz[®] Cover Glasses, 22 x 50 mm, precleaned: sc-24975
- PBS (Phosphate Buffered Saline), powder, 1 packet: sc-24947
- Formaldehyde, 37% formaldehyde solution, 25 ml: sc-203049
- Hydrogen Peroxide, 30% solution, 100 ml: sc-203336
- FCM Lysing solution: sc-3621
- FCM Fixation Buffer: sc-3622
- FCM Permeabilization Buffer: sc-3623
- FCM Wash Buffer: sc-3624
- Intracellular FCM System: sc-45063

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.