

v-Myb/c-Myb (SPM175): sc-56608

BACKGROUND

The Myb family of transcription factors regulates proliferation, differentiation and apoptosis of hematopoietic cells. The avian myeloblastosis viral (v-Myb) protein is nuclear and binds to specific DNA sequences. The gene encoding for v-Myb is oncogenic, causing monoblastic leukemia and transforming myelomonocytic cells by deregulating the expression of specific target genes. v-Myb functions as a transcriptional activator, and it can repress biologically relevant genes such as Ets-2, which promotes macrophage differentiation. The proto-oncogene c-Myb encodes for a nuclear protein that plays a role in transcriptional regulation and may be essential for hematopoietic cell proliferation. A-Myb is also member of the Myb family that is expressed in proliferating B cell centroblasts. Transgenic mice overexpressing A-Myb possess enhanced hyperplasia of the lymph nodes.

REFERENCES

1. Dvorakova, M., et al. 2001. An *ex vivo* model to study v-Myb-induced leukemogenicity. *Blood Cells Mol. Dis.* 27: 437-445.
2. Tahirov, T.H., et al. 2001. Crystals of ternary protein-DNA complexes composed of DNA-binding domains of c-Myb or v-Myb, C/EBP α or C/EBP β and tom-1A promoter fragment. *Acta Crystallogr. D, Biol. Crystallogr.* 57: 1655-1658.
3. Bartonek, P., et al. 2002. bFGF signaling and v-Myb cooperate in sustained growth of primitive progenitors. *Oncogene* 21: 400-410.
4. Wang, D.M. and Lipsick, J.S. 2002. Mutational analysis of the transcriptional activation domains of v-Myb. *Oncogene* 21: 1611-1615.
5. Bryja, V., et al. 2003. Lipoxigenase inhibitors v-Myb-transformed monoblasts BM2. *Prostaglandins Other Lipid Mediat.* 72: 131-145.
6. Nemajerova, A., et al. 2003. Trichostatin A suppresses transformation by the v-Myb oncogene in BM2 cells. *J. Hematother. Stem Cell Res.* 12: 225-235.
7. Chayka, O., et al. 2004. v-Myb mediates cooperation of a cell-specific enhancer with the mim-1 promoter. *Mol. Cell. Biol.* 25: 499-511.
8. Kanei-Ishii, C., et al. 2004. Differential sensitivity of v-Myb and c-Myb to Wnt-1-induced protein degradation. *J. Biol. Chem.* 279: 44582-44589.

CHROMOSOMAL LOCATION

Genetic locus: MYB (human) mapping to 6q23.3; Myb (mouse) mapping to 10 A3.

SOURCE

v-Myb/c-Myb (SPM175) is a mouse monoclonal antibody raised against recombinant full length v-Myb/c-Myb of human origin.

PRODUCT

Each vial contains 50 μ g IgG_{2a} in 0.5 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS


v-Myb/c-Myb (SPM175) is recommended for detection of both v-Myb and c-Myb of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Molecular Weight of c-Myb: 75 kDa.

Molecular Weight of v-Myb: 72 kDa.

Positive Controls: CCRF-CEM nuclear extract: sc-2146, Jurkat whole cell lysate: sc-2204 or Jurkat nuclear extract: sc-2132.

DATA



c-Myb (SPM175): sc-56608. Western blot analysis of c-Myb expression in Jurkat (A) and CCRF-CEM (B) nuclear extracts.

SELECT PRODUCT CITATIONS

1. Persson, M., et al. 2022. Rearrangements, expression, and clinical significance of MYB and MYBL1 in adenoid cystic carcinoma: a multi-institutional study. *Cancers* 14: 3691.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.



See **c-Myb (D-7): sc-74512** for c-Myb antibody conjugates, including AC, HRP, FITC, PE, and Alexa Fluor® 488, 546, 594, 647, 680 and 790.