

CD33 (CD33-4D3): sc-59020

BACKGROUND

CD33 is a type I transmembrane glycoprotein that is found on granulocyte and macrophage precursors in the bone marrow, but is absent from pluripotent stem cells. CD33 is also expressed on monocytes in peripheral blood. It is used as a marker to distinguish myelogenous leukemia cells from lymphoid or erythroid leukemias. CD33 may function as a sialic acid-dependent cell adhesion molecule.

REFERENCES

- Griffin, J.D., Linch, D., Sabbath, K., Larcom, P. and Schlossman, S.F. 1984. A monoclonal antibody reactive with normal and leukemic human myeloid progenitor cells. *Leuk. Res.* 8: 521-534.
- Favaloro, E.J., Bradstock, K.F., Kabral, A., Grimsley, P. and Berndt, M.C. 1987. Characterization of monoclonal antibodies to the human myeloid-differentiation antigen, "gp67" (CD-33). *Dis. Markers* 5: 215-225.
- Favaloro, E.J., Bradstock, K.F., Kabral, A., Grimsley, P., Zowtyj, H. and Zola, H. 1988. Further characterization of human myeloid antigens (gp160,95; gp150; gp67): investigation of epitopic heterogeneity and non-haemopoietic distribution using panels of monoclonal antibodies belonging to CD11b, CD13 and CD33. *Br. J. Haematol.* 69: 163-171.
- Knapp, W., Dorken, B., Gilks, W.R., Rieber, E.P., Schmidt, R.E., Stein, H. and Von dem Borne, A.E.G.K., eds. 1989. *Leucocyte Typing IV*. New York: Oxford University Press.
- Bradstock, K.F., Kirk, J., Grimsley, P.G., Kabral, A. and Hughes, W.G. 1989. Unusual immunophenotypes in acute leukaemias: incidence and clinical correlations. *Br. J. Haematol.* 72: 512-518.
- Robertson, M.J., Cunoosamy, M.P. and Clark, K.L. 1992. Effects of peptidase inhibition on angiotensin receptor agonist and antagonist potency in rabbit isolated thoracic aorta. *Br. J. Pharmacol.* 106: 166-172.
- Schlossman, S.F., Boumsell, L., Gilks, W., Harlan, J.M., Kishimoto, T., Morimoto, C., Ritz, J., Shaw, S., Silverstein, R., Springer, T., Tedder, T.F. and Todd, R.F., eds. 1995. *Leucocyte Typing V*. New York: Oxford University Press.
- Freeman, S.D., Kelm, S., Barber, E.K. and Crocker, P.R. 1995. Characterization of CD33 as a new member of the sialoadhesin family of cellular interaction molecules. *Blood* 85: 2005-2012.
- Kelm, S., Schauer, R. and Crocker, P.R. 1996. The sialoadhesins—a family of sialic-acid-dependent cellular recognition molecules within the immunoglobulin superfamily. *Glycoconj. J.* 13: 913-926.

CHROMOSOMAL LOCATION

Genetic locus: CD33 (human) mapping to 19q13.41.

SOURCE

CD33 (CD33-4D3) is a mouse monoclonal antibody raised against CD33 of human origin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Each vial contains 200 µg IgG_{2b} kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

CD33 (CD33-4D3) is recommended for detection of CD33 of human origin by flow cytometry (1 µg per 1 x 10⁶ cells).

Suitable for use as control antibody for CD33 siRNA (h): sc-42782, CD33 shRNA Plasmid (h): sc-42782-SH and CD33 shRNA (h) Lentiviral Particles: sc-42782-V.

Molecular Weight of CD33: 67 kDa.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.



See **CD33 (6C5/2): sc-53199** for CD33 antibody conjugates, including AC, HRP, FITC, PE, and Alexa Fluor[®] 488, 546, 594, 647, 680 and 790.